




Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.924 (0.858 to 0.990)	0.905 (0.839 to 0.971)
Intercept:	3.3 (0.5 to 6.0)	4.0 (1.3 to 6.7)
Std Err Est:	0.6	0.6

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9778	Points (Plotted/Total): 38/38
Bias: 0.1	Outliers: None
XMean ± SD: 41.1 ± 3.1	Degrees Freedom: 36
YMean ± SD: 41.2 ± 2.8	Scatter Plot Bounds: None
Std Dev Diff: 0.7	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	10 Feb 2003	10 Feb 2003
Rep SD:	1	1
Result Ranges:	36 to 48	36 to 47
Units:	g/L	g/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	41	41	0	S00014	42	42	0	S00027	40	41	1
S00002	37	38	1	S00015	44	44	0	S00028	40	40	0
S00003	44	43	-1	S00016	42	43	1	S00029	36	37	1
S00004	40	41	1	S00017	41	41	0	S00030	40	41	1
S00005	37	38	1	S00018	36	36	0	S00031	46	46	0
S00006	38	39	1	S00019	38	38	0	S00032	42	42	0
S00007	41	41	0	S00020	42	42	0	S00033	48	47	-1
S00008	43	44	1	S00021	42	41	-1	S00034	41	41	0
S00009	41	41	0	S00022	42	41	-1	S00035	43	43	0
S00010	36	36	0	S00023	41	40	-1	S00036	39	39	0
S00011	39	39	0	S00024	47	46	-1	S00037	42	42	0
S00012	37	38	1	S00025	38	38	0	S00038	46	47	1
S00013	43	43	0	S00026	45	45	0				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.