




## Alternate Method Comparison

X Method: BD

Y Method: Greiner



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	0.994 (0.978 to 1.011)	0.993 (0.977 to 1.010)
<b>Intercept:</b>	0.1 (-1.5 to 1.7)	0.2 (-1.4 to 1.9)
<b>Std Err Est:</b>	1.8	1.8

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R):	0.9987	Points (Plotted/Total):	40/40
Bias:	-0.4	Outliers:	None
XMean ± SD:	92.4 ± 34.8	Degrees Freedom:	38
YMean ± SD:	92.0 ± 34.6	Scatter Plot Bounds:	None
Std Dev Diff:	1.8		
SubRange Bounds:	None		

### Experiment Description

	X Method	Y Method
ExptDate:	10 Feb 2003	10 Feb 2003
Rep SD:	1	1
Result Ranges:	36 to 224	36 to 224
Units:	U/L	U/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_  
Signature

\_\_\_\_\_ Date

## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	84	84	0	S00015	126	126	0	S00029	60	60	0
S00002	86	85	-1	S00016	107	106	-1	S00030	96	92	-4
S00003	98	97	-1	S00017	109	115	6	S00031	155	153	-2
S00004	88	87	-1	S00018	62	61	-1	S00032	82	81	-1
S00005	111	111	0	S00019	106	107	1	S00033	36	36	0
S00006	81	81	0	S00020	65	65	0	S00034	102	102	0
S00007	64	64	0	S00021	131	133	2	S00035	85	82	-3
S00008	79	79	0	S00022	124	124	0	S00036	70	70	0
S00009	65	65	0	S00023	86	86	0	S00037	59	60	1
S00010	136	135	-1	S00024	68	67	-1	S00038	160	156	-4
S00011	56	56	0	S00025	77	77	0	S00039	66	68	2
S00012	97	97	0	S00026	103	98	-5	S00040	73	73	0
S00013	60	60	0	S00027	92	90	-2				
S00014	224	224	0	S00028	68	68	0				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.