


## Alternate Method Comparison


X Method: BD

Y Method: Greiner


**Scatter Plot**



**Bias**



**Percent Bias**



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	1.000 (0.995 to 1.004)	0.999 (0.995 to 1.004)
<b>Intercept:</b>	0.1 (-0.1 to 0.4)	0.1 (-0.1 to 0.4)
<b>Std Err Est:</b>	0.7	0.7

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R):	0.9999	Points (Plotted/Total):	40/40
Bias:	0.1	Outliers:	None
XMean ± SD:	30.4 ± 52.6	Degrees Freedom:	38
YMean ± SD:	30.6 ± 52.5	Scatter Plot Bounds:	None
Std Dev Diff:	0.7		
SubRange Bounds:	None		

### Experiment Description

	X Method	Y Method
ExptDate:	11 Feb 2003	11 Feb 2003
Rep SD:	1	1
Result Ranges:	9 to 348	9 to 348
Units:	U/L	U/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_

Signature

\_\_\_\_\_

Date

## Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	18	18	0	S00015	348	348	0	S00029	32	32	0
S00002	11	12	1	S00016	46	46	0	S00030	20	19	-1
S00003	28	28	0	S00017	13	14	1	S00031	14	13	-1
S00004	18	19	1	S00018	20	19	-1	S00032	22	22	0
S00005	18	18	0	S00019	14	15	1	S00033	10	10	0
S00006	42	43	1	S00020	49	49	0	S00034	9	9	0
S00007	16	16	0	S00021	19	19	0	S00035	31	30	-1
S00008	23	24	1	S00022	40	40	0	S00036	18	18	0
S00009	12	13	1	S00023	26	25	-1	S00037	38	39	1
S00010	23	23	0	S00024	14	14	0	S00038	13	12	-1
S00011	13	14	1	S00025	20	21	1	S00039	20	21	1
S00012	13	13	0	S00026	27	26	-1	S00040	34	35	1
S00013	11	12	1	S00027	20	19	-1				
S00014	16	16	0	S00028	38	38	0				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.