


Alternate Method Comparison

X Method: BD

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	0.996 (0.984 to 1.008)	0.995 (0.983 to 1.008)
Intercept:	0.3 (-0.2 to 0.7)	0.3 (-0.2 to 0.7)
Std Err Est:	1.0	1.0

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9993	Points (Plotted/Total):	40/40
Bias:	0.2	Outliers:	None
XMean ± SD:	25.6 ± 26.2	Degrees Freedom:	38
YMean ± SD:	25.8 ± 26.1	Scatter Plot Bounds:	None
Std Dev Diff:	1.0		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	10 Feb 2003	10 Feb 2003
Rep SD:	1	1
Result Ranges:	10 to 174	10 to 174
Units:	U/L	U/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	18	18	0	S00015	174	174	0	S00029	28	28	0
S00002	10	10	0	S00016	26	30	4	S00030	16	17	1
S00003	24	25	1	S00017	22	23	1	S00031	20	19	-1
S00004	12	11	-1	S00018	16	17	1	S00032	33	32	-1
S00005	26	25	-1	S00019	22	22	0	S00033	16	16	0
S00006	22	23	1	S00020	75	74	-1	S00034	14	14	0
S00007	17	18	1	S00021	21	21	0	S00035	21	22	1
S00008	21	21	0	S00022	28	29	1	S00036	18	19	1
S00009	14	15	1	S00023	21	21	0	S00037	24	24	0
S00010	36	35	-1	S00024	19	18	-1	S00038	14	15	1
S00011	17	16	-1	S00025	18	18	0	S00039	11	12	1
S00012	19	19	0	S00026	28	27	-1	S00040	21	22	1
S00013	17	17	0	S00027	16	16	0				
S00014	23	23	0	S00028	28	27	-1				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.