




Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.920 (0.837 to 1.003)	0.888 (0.805 to 0.971)
Intercept:	0.211 (0.018 to 0.404)	0.285 (0.094 to 0.477)
Std Err Est:	0.031	0.031

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9621	Points (Plotted/Total):	40/40
Bias:	0.026	Outliers:	None
XMean ± SD:	2.315 ± 0.120	Degrees Freedom:	38
YMean ± SD:	2.341 ± 0.111	Scatter Plot Bounds:	None
Std Dev Diff:	0.033		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	10 Feb 2003	10 Feb 2003
Rep SD:	1	1
Result Ranges:	2.05 to 2.58	2.11 to 2.62
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____
Signature

_____ Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	2.28	2.32	0.04	S00015	2.36	2.42	0.06	S00029	2.25	2.30	0.05
S00002	2.17	2.19	0.02	S00016	2.32	2.38	0.06	S00030	2.27	2.26	-0.01
S00003	2.23	2.29	0.06	S00017	2.58	2.62	0.04	S00031	2.40	2.44	0.04
S00004	2.09	2.12	0.03	S00018	2.36	2.32	-0.04	S00032	2.44	2.48	0.04
S00005	2.05	2.19	0.14	S00019	2.27	2.25	-0.02	S00033	2.31	2.36	0.05
S00006	2.26	2.27	0.01	S00020	2.21	2.27	0.06	S00034	2.38	2.41	0.03
S00007	2.33	2.36	0.03	S00021	2.18	2.20	0.02	S00035	2.47	2.49	0.02
S00008	2.29	2.31	0.02	S00022	2.15	2.22	0.07	S00036	2.48	2.50	0.02
S00009	2.33	2.31	-0.02	S00023	2.29	2.29	0.00	S00037	2.34	2.35	0.01
S00010	2.28	2.29	0.01	S00024	2.36	2.35	-0.01	S00038	2.50	2.48	-0.02
S00011	2.29	2.29	0.00	S00025	2.26	2.27	0.01	S00039	2.44	2.46	0.02
S00012	2.06	2.11	0.05	S00026	2.48	2.45	-0.03	S00040	2.45	2.48	0.03
S00013	2.44	2.47	0.03	S00027	2.25	2.31	0.06				
S00014	2.30	2.35	0.05	S00028	2.39	2.41	0.02				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.