




Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	0.963 (0.858 to 1.069)	0.912 (0.807 to 1.016)
Intercept:	1.9 (-9.3 to 13.1)	7.4 (-3.6 to 18.4)
Std Err Est:	1.1	1.1

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9442	Points (Plotted/Total):	40/40
Bias:	-1.9	Outliers:	None
XMean ± SD:	105.5 ± 3.4	Degrees Freedom:	38
YMean ± SD:	103.5 ± 3.2	Scatter Plot Bounds:	None
Std Dev Diff:	1.1		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	27 Jan 2003	27 Jan 2003
Rep SD:	1	1
Result Ranges:	100 to 112	96 to 110
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	103	102	-1	S00015	109	108	-1	S00029	106	104	-2
S00002	107	106	-1	S00016	108	105	-3	S00030	103	102	-1
S00003	107	106	-1	S00017	106	103	-3	S00031	105	104	-1
S00004	103	103	0	S00018	112	110	-2	S00032	112	110	-2
S00005	104	104	0	S00019	103	101	-2	S00033	101	101	0
S00006	108	105	-3	S00020	105	103	-2	S00034	105	103	-2
S00007	100	99	-1	S00021	105	103	-2	S00035	103	99	-4
S00008	103	101	-2	S00022	109	107	-2	S00036	107	106	-1
S00009	102	100	-2	S00023	101	100	-1	S00037	110	107	-3
S00010	102	101	-1	S00024	112	107	-5	S00038	104	103	-1
S00011	106	103	-3	S00025	104	101	-3	S00039	108	106	-2
S00012	111	109	-2	S00026	110	108	-2	S00040	102	100	-2
S00013	104	102	-2	S00027	103	100	-3				
S00014	100	96	-4	S00028	105	102	-3				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.