


Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.997 (0.864 to 1.130)	0.915 (0.784 to 1.045)
Intercept:	3.1 (0.1 to 6.1)	4.9 (2.0 to 7.8)
Std Err Est:	1.4	1.3

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9172	Points (Plotted/Total):	40/40
Bias:	3.0	Outliers:	None
XMean ± SD:	22.0 ± 3.3	Degrees Freedom:	38
YMean ± SD:	25.1 ± 3.3	Scatter Plot Bounds:	None
Std Dev Diff:	1.3		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	06 Feb 2003	06 Feb 2003
Rep SD:	1	1
Result Ranges:	15 to 30	18 to 33
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____
Signature

_____ Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	20	24	4	S00015	22	26	4	S00029	19	23	4
S00002	23	26	3	S00016	17	23	6	S00030	26	27	1
S00003	22	24	2	S00017	19	22	3	S00031	20	21	1
S00004	26	30	4	S00018	15	19	4	S00032	17	18	1
S00005	21	24	3	S00019	22	25	3	S00033	23	23	0
S00006	21	25	4	S00020	26	28	2	S00034	23	24	1
S00007	25	28	3	S00021	19	23	4	S00035	23	27	4
S00008	23	26	3	S00022	19	21	2	S00036	24	26	2
S00009	21	21	0	S00023	21	24	3	S00037	23	25	2
S00010	22	26	4	S00024	18	23	5	S00038	30	33	3
S00011	18	23	5	S00025	25	28	3	S00039	25	28	3
S00012	19	22	3	S00026	28	31	3	S00040	26	30	4
S00013	26	30	4	S00027	23	27	4				
S00014	18	22	4	S00028	23	27	4				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.