


Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	1.036 (0.987 to 1.086)	1.025 (0.976 to 1.074)
Intercept:	-13.2 (-31.1 to 4.7)	-9.5 (-27.3 to 8.4)
Std Err Est:	24.6	24.5

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9894	Points (Plotted/Total): 40/40
Bias: -1.4	Outliers: None
XMean ± SD: 325.4 ± 161.1	Degrees Freedom: 38
YMean ± SD: 324.0 ± 166.8	Scatter Plot Bounds: None
Std Dev Diff: 24.5	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	07 Feb 2003	07 Feb 2003
Rep SD:	1	1
Result Ranges:	159 to 915	109 to 964
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	314	323	9	S00015	261	208	-53	S00029	298	287	-11
S00002	246	257	11	S00016	503	511	8	S00030	223	200	-23
S00003	170	173	3	S00017	786	758	-28	S00031	266	263	-3
S00004	298	293	-5	S00018	252	244	-8	S00032	252	260	8
S00005	159	171	12	S00019	508	509	1	S00033	299	292	-7
S00006	299	266	-33	S00020	273	338	65	S00034	228	226	-2
S00007	915	964	49	S00021	232	218	-14	S00035	226	247	21
S00008	253	258	5	S00022	242	262	20	S00036	212	221	9
S00009	200	180	-20	S00023	293	315	22	S00037	256	270	14
S00010	653	645	-8	S00024	399	423	24	S00038	420	394	-26
S00011	379	359	-20	S00025	254	252	-2	S00039	270	290	20
S00012	491	494	3	S00026	267	248	-19	S00040	183	109	-74
S00013	264	275	11	S00027	261	256	-5				
S00014	470	481	11	S00028	242	221	-21				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.