


Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.989 (0.983 to 0.994)	0.989 (0.983 to 0.994)
Intercept:	-0.8 (-2.0 to 0.3)	-0.8 (-1.9 to 0.3)
Std Err Est:	2.4	2.4

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9999	Points (Plotted/Total):	40/40
Bias:	-2.6	Outliers:	None
XMean ± SD:	160.1 ± 151.0	Degrees Freedom:	38
YMean ± SD:	157.4 ± 149.3	Scatter Plot Bounds:	None
Std Dev Diff:	3.0		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	06 Feb 2003	06 Feb 2003
Rep SD:	1	1
Result Ranges:	55 to 785	55 to 777
Units:	umol/L	umol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	148	146	-2	S00015	152	152	0	S00029	97	97	0
S00002	122	120	-2	S00016	55	55	0	S00030	109	104	-5
S00003	132	129	-3	S00017	424	420	-4	S00031	517	508	-9
S00004	118	113	-5	S00018	208	208	0	S00032	545	537	-8
S00005	92	90	-2	S00019	72	74	2	S00033	97	101	4
S00006	74	74	0	S00020	102	97	-5	S00034	64	67	3
S00007	102	97	-5	S00021	79	79	0	S00035	108	104	-4
S00008	115	111	-4	S00022	296	293	-3	S00036	92	88	-4
S00009	141	139	-2	S00023	79	79	0	S00037	116	113	-3
S00010	69	69	0	S00024	115	109	-6	S00038	113	109	-4
S00011	125	120	-5	S00025	92	90	-2	S00039	88	86	-2
S00012	785	777	-8	S00026	222	219	-3	S00040	67	67	0
S00013	120	115	-5	S00027	85	83	-2				
S00014	157	150	-7	S00028	109	109	0				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.