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KOVA® System for Standardized Urinalysis

KOVA SYSTEM FOR STANDARDIZED URINALYSIS

Urinalysis, as currently performed in many laboratories, is carried out using various nonstandard procedures. These procedures vary from laboratory to laboratory and often the actual technique within the laboratory varies depending on the person performing the tests.

Sources of variation in conventional urinalysis:

- variable urine volumes
- different centrifugal conditions creating varying amounts of sediment for microscopic examination
- different amounts of sediment collected and suspended under the cover glass
- technique variations among individuals performing the procedure.

To standardize the urinalysis procedure, a constant specimen volume, centrifugal force and sediment volume must be maintained and a consistent method of microscopic examination and reporting of results should be used. Kova's KOVA System achieves this standardization by reducing variation, including technique differences among technicians.

INTENDED USE

The KOVA System offers a procedure and products that can be used to produce standardized results during routine urinalysis. Volume control, consistency and hygiene are provided from collection and transport to microscopic analysis of urine sediment. KOVA-Trol™ can be used for complete quality control of physical, chemical and microscopic examination test procedures.

ADVANTAGES

If the described procedure is followed consistently, one can use the values obtained in urinalysis with confidence. Clinicians can follow the progress and treatment of patients with certainty; any changes that occur outside the narrower limits that this system allows can be considered significant.

Laboratories may be compared and patients under observation can have their urinalysis done at different laboratories with comparable results.

KOVA SYSTEM AND SYSTEM COMPONENTS

Product Number	Product Description	Determinations Per Package
87153	KOVA System Super Pac 1000 w/Caps 100 KOVA Glasstic® Slide 10 (10 chambered), 1000 KOVA Petters, 1000 KOVA Super Tubes, 1000 KOVA Caps	1000
87154	KOVA System Super Pac 1000 100 KOVA Glasstic® Slide 10 (10 chambered), 1000 KOVA Petters, 1000 KOVA Super Tubes,	1000
87162	KOVA System Super Pac 1000 with Grids 100 KOVA Glasstic® Slide 10 (10 chambered) w/Grids, 1000 KOVA Petters, 1000 KOVA Super Tubes,	1000
87155	KOVA System Pac II 100 KOVA Slide II (4 chambered), 400 KOVA Petters, 400 KOVA Super Tubes,	400
87159	KOVA Mini Pack 5 KOVA Glasstic® Slide 10, 20 KOVA Petters, 20 KOVA Tube, 20 KOVA Caps	50
87156	KOVA System Value Pac 500 50 KOVA Glasstic Slide 10 with grids 500 KOVA Economy Tubes, 100 KOVA Caps	500
87158	KOVA System Value Pac 500 with Grids 50 KOVA Glasstic Slide 10 (10 chambered), 500 KOVA Petters, 500 Economy Tubes	500
87141	KOVA KO-LEC-PAC® 500 KOVA Super Tubes, 500 KOVA Caps, 500 KOVA Cups, 500 Labels and 5 Transport Racks	500
87100	KOVA Slide II with Grid for quantitation; 100 x 4 well slides; with each 1mm x 1mm grid square	400
87118	KOVA Slide II (without grid) 100 x 4 well slides	400
87146	KOVA Glasstic® Slide I0 100 x 10 well slides in crystal clear Acrylic	1000
87157	KOVA Glasstic® Slide I0 50 x 10 well slides in crystal clear Acrylic	500

87144	KOVA Glasstic® Slide I0 with Grid 100 x 10 well slides in crystal clear Plexiglas* with quantitation grids; each chamber contains 6.6 µl and has a 3 mm x 3 mm grid with fine divisions of 0.33 mm x 0.33 mm. The test procedure includes a method for quantitating cells per µl of patient samples.	1000
87137	KOVA Super Tube Graduated non-sterile disposable collection and centrifuge tubes made of high impact, unbreakable plastic to eliminate cracking or breaking during centrifugation.	500
87138	KOVA Economy Tube As above but in economical, break-resistant styrene plastic.	500
87135	KOVA Petter Disposable plastic transfer pipette designed to retain 1.0 ml of urine after centrifugation. The unique lock tip provides a one-step contamination-free decanting method.	500
87139	KOVA Cap Recommended to prevent spillage during transportation, as well as aerosol contamination during centrifugation.	500
87136	KOVA Decanting Rack Rack for decanting up to 10 specimens.	1 Rack
87116	KOVA Stain Modified Sternheimer Malbin stain, facilitates differentiation of formed elements.	3 x 25 mls

KOVA SYSTEM FOR QUALITY CONTROL - KOVA Controls

Human urine controls for physical, chemical and microscopic quality control.

Product Number	Description KOVA-Trol	Packaging
87329	KOVA-Trol I - High Abnormal Without Abnormal Urobilinogen Value Assignment	4 x 15 mls
87325	KOVA-Trol I - High Abnormal Without Abnormal Urobilinogen Value Assignment	4 x 60 mls
87326	KOVA-Trol I - High Abnormal Without Abnormal Urobilinogen Value Assignment	10 x 60 mls
87332	KOVA-Trol I - High Abnormal With Urobilinogen Value Assignment	4 x 60 mls
87333	KOVA-Trol I - High Abnormal With Urobilinogen Value Assignment	10 x 60 mls
87334	KOVA-Trol I - High Abnormal With Urobilinogen Value Assignment	4 x 15 mls
87130	KOVA-Trol II - Low Abnormal	4 x 15 mls
87128	KOVA-Trol II - Low Abnormal	10 x 60 mls
87331	KOVA-Trol III - Normal with hCG	4 x 15 mls
87327	KOVA-Trol III - Normal with hCG	4 x 60 mls
87328	KOVA-Trol III - Normal with hCG	10 x 60 mls

Simulated urine controls for physical, chemical and microscopic quality control.

Product Number	Description KOVA Liqua-Trol with Microscopics	Packaging
87122	KOVA Liqua-Trol Level II (Normal with hCG)	2 x 120ml
87176	KOVA Liqua-Trol Level I (Abnormal)	2 x 120ml
87123	KOVA Liqua-Trol Level II (Normal with hCG)	4 x 120ml
87177	KOVA Liqua-Trol Level I (Abnormal)	4 x 120ml

Product Number	Description KOVA Liqua-Trol without Microscopics	Packaging
87112	KOVA Liqua-Trol Level II (Normal with hCG) Control KOVA Liqua-Trol Level I (Abnormal) Control	3 x 15ml 3 x 15ml

SPECIMEN COLLECTION AND TRANSPORT

The KOVA System KO-LEC-PAC® is recommended for use in the following manner:

1. Label the KOVA Tube and give the patient a 3 ½ oz. KOVA Cup.
2. Instruct the patient to collect the voided urine in the KOVA Cup.
3. Transfer the urine specimen from the KOVA Cup to the KOVA Tube, filling it to the 12 ml graduation.
4. Secure the KOVA Cap on the KOVA Tube and place it in the KOVA Transport Rack for transportation and storage.
5. Deliver to the laboratory as soon as possible, preferably within two hours, but no more than four hours following specimen collection.


KOVA SYSTEM TEST PROCEDURE

- Check the specific gravity by placing one or two drops of urine in a temperature-compensated refractometer, or use a chemistry test strip containing a specific gravity parameter and record the results.
- Using reagent test strips, perform chemical testing according to the manufacturer's instructions. Record the observed results. KOVA-Trols I, II and III should be included in each batch to assure proper quality control of physical, chemical and microscopic test procedures.
- Centrifuge the KOVA Tubes (each containing 12ml of urine specimen or KOVA-Trol) at a relative centrifugal force (rcf) of 400 for five minutes; approximately 1500 revolutions per minute (rpm) with a 6-inch radius rotor. Formula used:

$$rcf = 28.38 (R) \left(\frac{N}{1000} \right)^2 R = \text{radius of rotor in inches}$$

N = revolutions per minute

The rotating radius is the distance measured from the rotor axis to the tip of the liquid inside the tubes at the greatest horizontal distance from the rotor axis.

- Remove the KOVA Tubes from the centrifuge being careful not to disturb or dislodge the sediment.
- Insert a KOVA Petter into the KOVA Tube. Push the KOVA Petter to the bottom of the KOVA Tube until it seats firmly (at the 1ml graduation).
- Decant and discard 11ml from the KOVA Tube while the KOVA Petter is locked in position in the KOVA Tube. This will retain 1ml of urine sediment at the bottom of the KOVA Tube.
- Withdraw the KOVA Petter from the KOVA Tube.
- Add one drop of KOVA Stain to the 1ml of urine sediment.
- Using the KOVA Petter, gently resuspend the sediment and stain until a homogeneous mixture is obtained.
- Withdraw a small sample of the urine sediment stain mixture by squeezing the bulb of the KOVA Petter.
- Transfer the sediment mixture to the KOVA Slide by placing one drop in the cut-out notch of each chamber. When chambers 1-5 are on the top row, the notch is at the top left corner of chambers, when chambers 6-10 are on the top row, the notch is at the top right corner of chambers. The chamber will fill by capillary action. Avoid touching the V-shaped barrier between the chambers while dispensing fluid. Incorrect positioning in dispensing may cause overflowing from one chamber to the next.
- Remove any excess specimen remaining on the open recessed area by touching the open edge with absorbent material.
- Place the KOVA Slide on a microscopic stage under the objective lens.
- Scan the slide chamber under low power magnification (10X eyepiece/10X objective) to enumerate casts. Enumerate all other formed elements under high power magnification (10X eyepiece/40X objective). Do not reuse KOVA products.

EXPECTED VALUES - MICROSCOPY†

Analyte	Normal	Abnormal	Reporting Results
WBC	0-5/HPF	> 5/HPF	Numbers/HPF
RBC	0-3/HPF	> 3/HPF	Numbers/HPF
Epithelial Cells	0	Any (other than squamous)	Numbers/HPF
Crystals	0-3/HPF (non-pathogenic)	> 3 Any abnormal	Numbers/HPF
Yeasts	0	Any	1 + to 4 +/HPF
Trichomonads	0	Any	1 + to 4 +/HPF
Casts	0	Any especially > 1 hyaline cast/LPF	Numbers/LPF
Bacteria	0-5/HPF	> 5/HPF	1 + to 4 +/HPF
Fat	0	Oval fat bodies or free fat	1 + to 4 +/HPF

† Bernard Statland, MLO. p 13-14; Jan. 1985

References for General Information

- Bradley, G.M., Benson, E.S., Todd-Sanford Clinical Diagnosis by Laboratory Methods, 15th Edition, Phila. Saunders, 1974.
- Kurtzman, N.A. and Rogers, P.W. (1974). A Handbook of Urinalysis and Urinary Sediment. Chas. C Thomas, Springfield, IL.
- Little, P.J. (1962). Urinary white-cell excretion. Lancet. pp. 1149-1151.
- Little, P.F. (1964). A comparison of the urinary white cell concentrations with the white cell excretion rate. Brit. J. Urol. 36, 360-363.
- Thomas, M.(1971). A rapid slide method of urine cell counts. Med. Lab Technol. 28, 38-39.
- Moore, T., Hira, N.R., and Stirland, R.M. (1965). Differential urethrovesical urinary cell count. Lancet. pp. 626-627.
- Siegle, M.D., Lab Med., 12:781, 1981.
- Sternheimer, R. and Malbin, B. (1951). The clinical recognition of pyelonephritis with a new stain for urinary sediments. Am. J. of Med., 11:312-323.
- Muschetta, P.A. and Waters , Jr. F.O. (1962). Manual of Medical Laboratory Techniques. Herbert-Spence, Inc. New York, N.Y., Second Edition, pp 44-45.
- Lippman, R. W. (1957). Urine and the Urinary Sediment. Chas. C Thomas, Springfield, IL.
- Dudas, H.C., Lab Med. 12:765. 1981.
- Weller, J.M. and Greene, J.A. (1966). Examination of the Urine. Meredith Publishing Co., New York.

KOVA, KO-LEC-PAC and Glasstic are registered trademarks of Kova International, Inc., Garden Grove, CA, USA.

KOVA-Trol is a trademark of Kova International, Inc., Garden Grove, CA, USA.

U.S. PATENT NUMBERS

4,022,576
4,563,332
4,299,441
D258,388
4,637,693
4,937,415
4,997,266
RE 33,826

KOVA
INTERNATIONAL
www.kovaintl.com

Kova International, Inc.
7272 Chapman Avenue, Suite B
Garden Grove, California 92841
UNITED STATES
Tel: +1 855 217 6399
Fax: +1 714 908 7945

Advena Ltd.
Tower Business Centre
2nd Flr, Tower Street
Swatar, Bkr 4013
Malta

Symbols and Translations

Symbol	English	French	German	Italian	Spanish	Greek	Swedish
	Batch / Lot code	Code du lot	Chargenbezeichnung	Numero di lotto	Código de lote	Παρτίδα / Κωδικός Παραγωγής	Batch-/lolkod
	Expiry / Use by	Date d'expiration	Verwendbar bis	Utilizzare entro	Fecha de caducidad	Λήξη/ Χρήση έως	Utgångsdatum/förbrukas senast
	Manufactured by	Fabriqué par	Hergestellt von	Fabbricato da	Fabricado por	Κατασκευασμένο από	Tillverkad av
	Store at	Conserver à	Lagerung bei	Conservare a	Conserver a	Αποθηκεύσατε στους	Förvaras vid
	See instructions for use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Συμβουλευτείτε τις οδηγίες χρήσης	Se bruksanvisningen
	For in vitro diagnostic use only	Pour utilisation diagnostic in vitro uniquement	In-Vitro-Diagnoslikum	Per uso diagnostico in vitro	Sólo para el uso diagnóstico in vitro	Προορισμένο για διάγνωση In Vitro	Endast för in vitro-diagnostisk användning
	Negative control	Contrôle négatif	Negativ Kontrolle	Controllo negativo	Control negativo	Αρνητικός ορός ελέγχου	Negativ kontroll
	Positive control	Contrôle positif	Positiv Kontrolle	Controllo positivo	Control positivo	Θετικός ορός ελέγχου	Positiv kontroll
	Catalogue number	numéro de produit	Bestellnummer	Numero di catalogo	Número de catálogo	Αριθμός Καταλόγου	Katalognummer
	Authorized representative	Représentant autorisé	Vertrieb	Distributore Autorizzato	Representante Autorizado	Εξουσιοδοτημένος αντιπρόσωπος	Auktoriserat ombud
	Contains sufficient for <n> tests	Contenu suffisant pour <n> tests	Enthält ausreichend Reagenzien für <n> Tests	Contenido suficiente per <n> test	Conviene lo necesario para <n> test	Περιέχει επαρκή ποσότητα για <n> ελέγχους	Räcker till <n> tester
	Do not reuse	Ne pas réutiliser.	Nicht wiederverwenden.	Non riutilizzare.	No reutilizar.	Μην επαναχρησιμοποιείτε	Återanvänds ej
	Warning	Avertissement	Warnung	Avvertimento	Advertencia	προειδοποίηση	Varning

P/N 97120-9