Operation Manual

Greiner Bio-One Microplate Centrifuge
REF 846050 (US) / 846060 (UK) / 846070 (EU)

Introduction:

Thank you for purchasing the Greiner Bio-One Microplate Centrifuge. Please read this manual thoroughly prior to operating the instrument. Using the device inconsistently with the instructions may result in damaging the device.

Text labeled with this sign contains safety notes. Read these safety precautions before using the device.

Intended Use:

The Microplate Centrifuge is intended only for indoor use and for separating aqueous solutions and suspensions of various densities predominantly in research laboratories. It is designed for quick spin downs of approved PCR plates and microtiter plates.

Specifications:

Speed:	2550 rpm / 600 x g
Capacity:	2 x PCR Plates (or microplates)
Dimensions (Width/Depth/Height):	9.2 (W) x 10.2 (D) x 7.75 (H) in. 23 (W) x 26 (D) x 19.7 (H) cm
Electrical:	120 V/60 Hz, 230 V/50 Hz
Max. Radius:	3.15 in. / 8.0 cm

Operating Environment:	4 °C to 45 °C 39.2 °F to 113 °F
Weight:	9 lbs. / 4,1 kg

Installation:

Place the Microplate Centrifuge on a stable, horizontal and level surface nearby an electrical outlet.

In accordance with the safety regulations of IEC 1010-2-020 a safety distance of 30 cm should be observed around the centrifuge during operation. The Centrifuge should not be near a heat or water source or in direct sunlight. Before plugging in the centrifuge, compare your power supply with the electrical requirements. Before starting the centrifuge for the first time, make sure the rotor is securely attached.

Safety Precautions:

The rotor must always be securely attached. If the centrifuge makes unusual noises when started, switch off the device immediately.

The rotor must only be loaded symmetrically (see also the section Loading the Rotor). Do not move the centrifuge during the run.

Do not use centrifuges that have not been correctly installed or repaired (see also Service and Contact)

The Microplate Centrifuge should only be used for the specified applications. It must not be operated in hazardous or flammable environments and must not be used to centrifuge explosive or highly reactive substances.

When handling toxic, radioactive or pathogenic liquids, observe the safety regulations of the country in question. If such liquids are spilled, the centrifuge must be cleaned carefully and properly.

Before using cleaning or decontamination methods others than those indicated in the manual, contact the manufacturer to ensure that the intended methods will not damage the centrifuge.

Plates and PCR strips should be visually inspected for material damage prior to centrifugation. Damaged products must not be centrifuged.

When moving the centrifuge from a cold room to a nor-

mal lab, allow it to warm up in the lab, but do not plug in the centrifuge in order to prevent damage caused by condensation.

Operation:

- Attach the power cord and place the power switch into the on "I" position.
- Open the lid and load your samples into the rotor. Always ensure a balanced load.

NOTE: The centrifuge should never be started with one plate or two plates that are not equally balanced. Additional information can be found in the following section of this manual Loading the Rotor.

- Close the lid. The rotor quickly accelerates to the maximum speed of 2550 rpm / 600 x g.
- When the desired run time is completed, push on the lid button to open the lid. The rotor quickly brakes to a stop and the samples can be retrieved.

NEVER attempt to remove samples until the rotor has come to a complete stop.

Loading the Rotor:

To ensure safe operation and long life of the instrument, samples must always be loaded into the rotor in a balanced fashion.

Plates must always be loaded in the center of the carrier. When loading non-skirted plates, it is important to line up the wells of the plate with those shown on the plate carrier. Starting the centrifuge with the plates not properly centered will result in an unbalanced load and increased vibration.

Example of properly balanced rotor:

Example of improperly balanced rotor:

Cleaning and Maintenance:

The outside of the Microplate Centrifuge can be cleaned with a moist cloth and a mild noncorrosive detergent. Do not immerse the centrifuge in liquid or pour liquids over it. Excessive liquids should be avoided and liquids should not come in contact with the motor. Make sure that the centrifuge is disconnected from the main supply during cleaning. The centrifuge rotor can be removed for cleaning by removing the center thumb screw and lifting the rotor upward. Once removed, the rotor can be cleaned with 70 % isopropanol/water mixture or sterilised in an autoclave at 121 °C / 250 °F for 20 minutes.

After cleaning, ensure that all parts are dried thoroughly before attempting to operate the unit. Allow the rotor once to accelerate to the maximum speed.

Should a spill of toxic, radioactive or pathogenic liquids occur within the rotor or the rotor chamber, the centrifuge must be decontaminated thoroughly. Do not use any cleaning or decontamination methods except those recommended in this manual.

Service and Contact:

This product is warranted to be free of defects in material and workmanship for a period of 2 years from the date of purchase. The warranty is valid only if the product is used in its intended purpose and within the guidelines specified in this instruction manual. In the event that service or technical support is required, please contact your nearest Greiner Bio-One office. When returning the centrifuge, ensure that the device is fully decontaminated and does not present any kind of health risk to our service personnel.

Disposal:

In the case the product is to be disposed of, the relevant legal regulations are to be observed. The disposal of electrical devices is regulated within the European Community by national regulations based on ED directive 2002/96/EC on waste electrical and electronic equipment (WEEE). As disposal regulations within EU may vary from country to country, please contact your supplier if necessary.

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Greiner Bio-One Co. Ltd.
Greiner Bio-One GmbH
Greiner Bio-One B.V.
Greiner Bio-One Ltd.
Greiner Bio-One North America Inc.