


Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	1.115 (1.045 to 1.185)	1.092 (1.022 to 1.162)
Intercept:	-15.3 (-46.1 to 15.4)	-7.9 (-38.5 to 22.6)
Std Err Est:	64.2	63.8

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9815	Points (Plotted/Total): 40/40
Bias: 22.0	Outliers: None
XMean ± SD: 325.2 ± 295.6	Degrees Freedom: 38
YMean ± SD: 347.2 ± 329.0	Scatter Plot Bounds: None
Std Dev Diff: 68.6	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	19 Dec 2002	19 Dec 2002
Rep SD:	1	1
Result Ranges:	67 to 1282	20 to 1397
Units:	pmol/L	pmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____
Signature

_____ Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	923	1126	203	S00015	375	302	-73	S00029	145	142	-3
S00002	984	1036	52	S00016	858	771	-87	S00030	299	401	102
S00003	204	196	-8	S00017	160	158	-2	S00031	284	155	-129
S00004	520	575	55	S00018	446	417	-29	S00032	269	267	-2
S00005	967	1109	142	S00019	115	216	101	S00033	319	350	31
S00006	666	816	150	S00020	163	171	8	S00034	97	154	57
S00007	1282	1397	115	S00021	109	127	18	S00035	115	98	-17
S00008	283	341	58	S00022	175	106	-69	S00036	80	156	76
S00009	145	237	92	S00023	214	190	-24	S00037	101	20	-81
S00010	297	368	71	S00024	229	190	-39	S00038	89	121	32
S00011	141	190	49	S00025	478	501	23	S00039	67	85	18
S00012	197	166	-31	S00026	104	107	3	S00040	89	62	-27
S00013	338	346	8	S00027	242	226	-16				
S00014	116	156	40	S00028	323	337	14				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.