




Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	1.016 (0.984 to 1.048)	1.011 (0.980 to 1.043)
Intercept:	-0.29 (-1.36 to 0.78)	-0.20 (-1.27 to 0.87)
Std Err Est:	2.72	2.72

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9955	Points (Plotted/Total): 40/40
Bias: 0.02	Outliers: None
XMean ± SD: 19.77 ± 27.93	Degrees Freedom: 38
YMean ± SD: 19.80 ± 28.38	Scatter Plot Bounds: None
Std Dev Diff: 2.70	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	03 Feb 2003	03 Feb 2003
Rep SD:	1	1
Result Ranges:	1.0 to 106.0	0.7 to 113.5
Units:	pmol/L	pmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	58.3	57.7	-0.6	S00015	72.2	65.4	-6.8	S00029	2.5	2.4	-0.1
S00002	4.7	4.8	0.1	S00016	3.2	2.9	-0.3	S00030	7.0	6.4	-0.6
S00003	4.4	4.6	0.2	S00017	1.5	1.5	0.0	S00031	33.5	35.9	2.4
S00004	55.2	51.4	-3.8	S00018	51.9	49.5	-2.4	S00032	1.9	1.6	-0.3
S00005	3.5	3.4	-0.1	S00019	1.7	1.8	0.1	S00033	2.8	2.6	-0.2
S00006	1.7	1.3	-0.4	S00020	2.6	2.6	0.0	S00034	61.1	57.9	-3.2
S00007	2.1	1.6	-0.5	S00021	1.0	0.7	-0.3	S00035	17.3	21.2	3.9
S00008	6.7	7.0	0.3	S00022	4.5	4.5	0.0	S00036	5.4	5.0	-0.4
S00009	6.5	5.7	-0.8	S00023	49.0	41.2	-7.8	S00037	4.7	6.1	1.4
S00010	1.6	1.4	-0.2	S00024	76.3	82.9	6.6	S00038	1.2	1.5	0.3
S00011	2.0	1.8	-0.2	S00025	4.7	4.3	-0.4	S00039	2.8	2.8	0.0
S00012	58.7	61.6	2.9	S00026	106.0	113.5	7.5	S00040	2.1	2.2	0.1
S00013	1.9	1.6	-0.3	S00027	57.0	60.3	3.3				
S00014	5.8	6.0	0.2	S00028	4.0	5.3	1.3				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.