


Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.968 (0.964 to 0.972)	0.968 (0.964 to 0.972)
Intercept:	0.3 (-0.4 to 0.9)	0.3 (-0.4 to 0.9)
Std Err Est:	2.0	2.0

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9999	Points (Plotted/Total): 40/40
Bias: -1.6	Outliers: Not tested
XMean ± SD: 59.7 ± 154.8	Degrees Freedom: 38
YMean ± SD: 58.1 ± 149.8	Scatter Plot Bounds: None
Std Dev Diff: 5.3	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	08 Jan 2003	08 Jan 2003
Rep SD:	1	1
Result Ranges:	7 to 893	7 to 869
Units:	U/L	U/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	24	23	-1	S00015	40	40	0	S00029	26	26	0
S00002	18	18	0	S00016	266	260	-6	S00030	26	25	-1
S00003	31	29	-2	S00017	18	18	0	S00031	22	21	-1
S00004	23	23	0	S00018	12	13	1	S00032	19	18	-1
S00005	30	30	0	S00019	18	18	0	S00033	54	54	0
S00006	15	15	0	S00020	19	19	0	S00034	10	10	0
S00007	26	25	-1	S00021	14	13	-1	S00035	21	20	-1
S00008	15	15	0	S00022	15	15	0	S00036	11	10	-1
S00009	32	32	0	S00023	24	25	1	S00037	19	18	-1
S00010	7	7	0	S00024	13	13	0	S00038	21	21	0
S00011	893	869	-24	S00025	13	14	1	S00039	13	13	0
S00012	26	25	-1	S00026	23	22	-1	S00040	12	12	0
S00013	24	23	-1	S00027	39	38	-1				
S00014	434	410	-24	S00028	23	23	0				

Values with an "X" were excluded from the calculations.