


Alternate Method Comparison

X Method: **BD**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	1.017 (0.986 to 1.048)	1.013 (0.982 to 1.044)
Intercept:	-0.09 (-0.28 to 0.10)	-0.06 (-0.25 to 0.13)
Std Err Est:	0.23	0.23

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R):	0.9956	Points (Plotted/Total):	41/41
Bias:	0.01	Outliers:	None
XMean ± SD:	5.68 ± 2.32	Degrees Freedom:	39
YMean ± SD:	5.69 ± 2.36	Scatter Plot Bounds:	None
Std Dev Diff:	0.22		
SubRange Bounds:	None		

Experiment Description

	X Method	Y Method
ExptDate:	27 Jan 2003	27 Jan 2003
Rep SD:	1	1
Result Ranges:	3.5 to 18.2	3.6 to 18.2
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____
Signature

_____ Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	5.7	5.6	-0.1	S00015	5.8	5.7	-0.1	S00029	7.8	7.8	0.0
S00002	5.1	5.0	-0.1	S00016	4.7	4.7	0.0	S00030	4.6	4.7	0.1
S00003	3.9	3.9	0.0	S00017	4.5	4.4	-0.1	S00031	4.7	4.6	-0.1
S00004	4.7	4.9	0.2	S00018	5.5	5.3	-0.2	S00032	4.6	4.5	-0.1
S00005	18.2	18.2	0.0	S00019	6.5	6.6	0.1	S00033	5.1	4.9	-0.2
S00006	5.8	5.8	0.0	S00020	5.0	5.1	0.1	S00034	5.5	5.5	0.0
S00007	5.1	5.2	0.1	S00021	4.4	4.4	0.0	S00035	6.5	6.6	0.1
S00008	7.5	7.5	0.0	S00022	5.2	5.4	0.2	S00036	5.6	5.7	0.1
S00009	3.5	3.6	0.1	S00023	8.0	8.1	0.1	S00037	5.5	5.5	0.0
S00010	5.5	5.4	-0.1	S00024	4.6	4.7	0.1	S00038	4.8	4.9	0.1
S00011	5.3	5.2	-0.1	S00025	5.5	5.4	-0.1	S00039	4.5	4.5	0.0
S00012	4.7	4.7	0.0	S00026	9.4	9.6	0.2	S00040	4.7	3.8	-0.9
S00013	4.4	4.4	0.0	S00027	4.5	4.6	0.1	S00041	7.0	7.9	0.9
S00014	4.3	4.3	0.0	S00028	4.7	4.7	0.0				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.