


## Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	0.984 (0.841 to 1.126)	0.890 (0.751 to 1.029)
<b>Intercept:</b>	-0.01 (-0.63 to 0.61)	0.40 (-0.21 to 1.00)
<b>Std Err Est:</b>	0.23	0.22

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R):	0.9028	Points (Plotted/Total):	40/40
Bias:	-0.08	Outliers:	None
XMean ± SD:	4.31 ± 0.52	Degrees Freedom:	38
YMean ± SD:	4.23 ± 0.51	Scatter Plot Bounds:	None
Std Dev Diff:	0.23		
SubRange Bounds:	None		

### Experiment Description

	X Method	Y Method
ExptDate:	03 Feb 2003	03 Feb 2003
Rep SD:	1	1
Result Ranges:	3.5 to 5.6	3.4 to 5.6
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_

Signature

\_\_\_\_\_

Date

## Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	3.8	3.7	-0.1	S00015	3.6	3.6	0.0	S00029	4.5	4.4	-0.1
S00002	4.3	4.2	-0.1	S00016	3.8	3.8	0.0	S00030	4.3	4.1	-0.2
S00003	4.5	4.3	-0.2	S00017	4.0	4.1	0.1	S00031	3.5	3.5	0.0
S00004	4.5	4.3	-0.2	S00018	5.6	5.6	0.0	S00032	5.4	5.3	-0.1
S00005	4.2	4.2	0.0	S00019	4.2	4.3	0.1	S00033	4.8	4.7	-0.1
S00006	4.8	4.4	-0.4	S00020	5.4	4.4	-1.0	S00034	4.5	4.4	-0.1
S00007	3.9	3.8	-0.1	S00021	3.7	3.8	0.1	S00035	3.6	3.4	-0.2
S00008	3.7	3.6	-0.1	S00022	4.7	4.6	-0.1	S00036	4.2	4.1	-0.1
S00009	4.3	4.5	0.2	S00023	4.4	4.4	0.0	S00037	4.5	4.4	-0.1
S00010	3.6	3.4	-0.2	S00024	4.0	4.0	0.0	S00038	4.0	4.0	0.0
S00011	4.1	3.7	-0.4	S00025	4.4	4.4	0.0	S00039	4.5	4.2	-0.3
S00012	4.8	5.1	0.3	S00026	3.9	4.4	0.5	S00040	4.8	5.0	0.2
S00013	4.0	3.9	-0.1	S00027	4.7	4.7	0.0				
S00014	4.9	4.8	-0.1	S00028	3.8	3.6	-0.2				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.