




Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	0.996 (0.953 to 1.040)	0.988 (0.945 to 1.031)
Intercept:	-11.1 (-20.2 to -2.1)	-9.4 (-18.4 to -0.4)
Std Err Est:	8.6	8.6

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9913	Points (Plotted/Total): 40/40
Bias: -11.8	Outliers: None
XMean ± SD: 199.3 ± 64.5	Degrees Freedom: 38
YMean ± SD: 187.5 ± 64.2	Scatter Plot Bounds: None
Std Dev Diff: 8.5	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	08 Jan 2003	08 Jan 2003
Rep SD:	1	1
Result Ranges:	108 to 364	106 to 342
Units:	U/L	U/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____
Signature

_____ Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	213	194	-19	S00015	232	225	-7	S00029	175	159	-16
S00002	169	164	-5	S00016	134	119	-15	S00030	155	133	-22
S00003	229	207	-22	S00017	326	322	-4	S00031	203	182	-21
S00004	160	161	1	S00018	162	151	-11	S00032	206	202	-4
S00005	166	156	-10	S00019	108	106	-2	S00033	364	342	-22
S00006	145	130	-15	S00020	162	156	-6	S00034	207	191	-16
S00007	174	160	-14	S00021	147	143	-4	S00035	228	199	-29
S00008	202	189	-13	S00022	190	186	-4	S00036	150	139	-11
S00009	328	342	14	S00023	155	148	-7	S00037	229	215	-14
S00010	154	138	-16	S00024	214	190	-24	S00038	131	128	-3
S00011	167	154	-13	S00025	182	172	-10	S00039	159	133	-26
S00012	208	201	-7	S00026	362	336	-26	S00040	160	152	-8
S00013	156	142	-14	S00027	151	143	-8				
S00014	278	268	-10	S00028	332	322	-10				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.