


Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.976 (0.931 to 1.020)	0.967 (0.923 to 1.011)
Intercept:	0.027 (-0.009 to 0.062)	0.034 (-0.002 to 0.070)
Std Err Est:	0.018	0.018

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9905	Points (Plotted/Total): 40/40
Bias: 0.007	Outliers: None
XMean ± SD: 0.803 ± 0.129	Degrees Freedom: 38
YMean ± SD: 0.810 ± 0.126	Scatter Plot Bounds: None
Std Dev Diff: 0.018	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	07 Feb 2003	07 Feb 2003
Rep SD:	1	1
Result Ranges:	0.56 to 1.24	0.57 to 1.23
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____

Signature

Date

Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	0.72	0.72	0.00	S00015	0.62	0.63	0.01	S00029	0.87	0.86	-0.01
S00002	0.95	0.98	0.03	S00016	0.80	0.82	0.02	S00030	0.81	0.80	-0.01
S00003	0.97	0.95	-0.02	S00017	1.01	1.05	0.04	S00031	1.24	1.23	-0.01
S00004	0.90	0.93	0.03	S00018	0.79	0.80	0.01	S00032	0.97	0.96	-0.01
S00005	0.85	0.86	0.01	S00019	0.75	0.74	-0.01	S00033	0.89	0.91	0.02
S00006	0.80	0.79	-0.01	S00020	0.80	0.80	0.00	S00034	0.72	0.73	0.01
S00007	0.94	0.90	-0.04	S00021	0.64	0.66	0.02	S00035	0.85	0.84	-0.01
S00008	0.71	0.72	0.01	S00022	0.81	0.82	0.01	S00036	0.70	0.71	0.01
S00009	0.64	0.69	0.05	S00023	0.70	0.71	0.01	S00037	0.67	0.68	0.01
S00010	0.64	0.65	0.01	S00024	0.74	0.76	0.02	S00038	0.86	0.88	0.02
S00011	0.56	0.57	0.01	S00025	0.92	0.94	0.02	S00039	0.77	0.75	-0.02
S00012	0.91	0.89	-0.02	S00026	0.80	0.82	0.02	S00040	0.84	0.84	0.00
S00013	0.79	0.82	0.03	S00027	0.74	0.74	0.00				
S00014	0.67	0.68	0.01	S00028	0.77	0.78	0.01				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.