




Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner



Regression Analysis

	Deming	Regular
Slope:	0.991 (0.869 to 1.112)	0.922 (0.803 to 1.042)
Intercept:	-0.7 (-18.4 to 16.9)	9.2 (-8.1 to 26.6)
Std Err Est:	1.1	1.1

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9306	Points (Plotted/Total): 40/40
Bias: -2.1	Outliers: None
XMean ± SD: 145.6 ± 3.0	Degrees Freedom: 38
YMean ± SD: 143.5 ± 2.9	Scatter Plot Bounds: None
Std Dev Diff: 1.1	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	03 Feb 2003	03 Feb 2003
Rep SD:	1	1
Result Ranges:	138 to 154	134 to 152
Units:	mmol/L	mmol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____
Signature

_____ Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

Alternate Method Comparison

X Method: BD - SST

Y Method: Greiner

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	142	141	-1	S00015	147	146	-1	S00029	145	142	-3
S00002	148	146	-2	S00016	144	141	-3	S00030	146	144	-2
S00003	146	144	-2	S00017	146	144	-2	S00031	144	143	-1
S00004	147	147	0	S00018	146	144	-2	S00032	144	142	-2
S00005	144	143	-1	S00019	143	141	-2	S00033	142	142	0
S00006	146	143	-3	S00020	146	143	-3	S00034	143	142	-1
S00007	145	144	-1	S00021	143	141	-2	S00035	147	142	-5
S00008	146	143	-3	S00022	146	145	-1	S00036	148	146	-2
S00009	141	140	-1	S00023	142	142	0	S00037	150	147	-3
S00010	143	141	-2	S00024	151	147	-4	S00038	149	147	-2
S00011	146	143	-3	S00025	145	143	-2	S00039	148	146	-2
S00012	149	147	-2	S00026	154	152	-2	S00040	144	141	-3
S00013	147	145	-2	S00027	143	140	-3				
S00014	138	134	-4	S00028	148	145	-3				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.