




## Alternate Method Comparison

**X Method: BD - SST**

**Y Method: Greiner**



### Regression Analysis

	Deming	Regular
<b>Slope:</b>	1.021 (0.993 to 1.048)	1.017 (0.989 to 1.045)
<b>Intercept:</b>	0.0 (-0.3 to 0.3)	0.0 (-0.3 to 0.3)
<b>Std Err Est:</b>	0.4	0.4

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R): 0.9966	Points (Plotted/Total): 40/40
Bias: 0.2	Outliers: None
XMean ± SD: 9.9 ± 5.2	Degrees Freedom: 38
YMean ± SD: 10.1 ± 5.3	Scatter Plot Bounds: None
Std Dev Diff: 0.4	
SubRange Bounds: None	

### Experiment Description

	X Method	Y Method
ExptDate:	19 Dec 2002	19 Dec 2002
Rep SD:	1	1
Result Ranges:	3 to 25	3 to 26
Units:	ug/L	ug/L
Analyst:	DIS	DIS
Comment:		

Accepted by: \_\_\_\_\_

Signature

\_\_\_\_\_

Date

Prepared for: Biochemistry Division -- Ottawa Civic Hospital

By: Technical Marketing -- greiner-bio-one

## Alternate Method Comparison

**X Method: BD - SST**

**Y Method: Greiner**

### Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	11	11	0	S00015	5	5	0	S00029	11	12	1
S00002	12	13	1	S00016	17	16	-1	S00030	9	9	0
S00003	17	17	0	S00017	15	15	0	S00031	13	13	0
S00004	3	3	0	S00018	9	9	0	S00032	5	5	0
S00005	5	5	0	S00019	13	14	1	S00033	10	10	0
S00006	7	7	0	S00020	23	23	0	S00034	6	6	0
S00007	6	7	1	S00021	5	5	0	S00035	15	15	0
S00008	25	26	1	S00022	15	15	0	S00036	10	10	0
S00009	18	19	1	S00023	10	10	0	S00037	10	10	0
S00010	9	9	0	S00024	4	4	0	S00038	14	15	1
S00011	6	6	0	S00025	8	8	0	S00039	7	7	0
S00012	8	8	0	S00026	5	5	0	S00040	6	6	0
S00013	8	8	0	S00027	5	6	1				
S00014	6	6	0	S00028	5	5	0				

Values with an "X" were excluded from the calculations. Outliers "O" were also excluded.