


Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**



Regression Analysis

	Deming	Regular
Slope:	0.951 (0.942 to 0.960)	0.951 (0.942 to 0.960)
Intercept:	0.2 (-0.1 to 0.6)	0.2 (-0.1 to 0.6)
Std Err Est:	1.0	1.0

95% Confidence Intervals are shown in parentheses

Supporting Statistics

Corr Coef (R): 0.9996	Points (Plotted/Total): 40/40
Bias: -0.6	Outliers: Not tested
XMean ± SD: 16.9 ± 35.8	Degrees Freedom: 38
YMean ± SD: 16.3 ± 34.0	Scatter Plot Bounds: None
Std Dev Diff: 2.0	
SubRange Bounds: None	

Experiment Description

	X Method	Y Method
ExptDate:	08 Jan 2003	08 Jan 2003
Rep SD:	1	1
Result Ranges:	3 to 191	3 to 180
Units:	umol/L	umol/L
Analyst:	DIS	DIS
Comment:		

Accepted by: _____
Signature

_____ Date

Alternate Method Comparison

X Method: **BD - SST**

Y Method: **Greiner**

Experimental Results

Specimen	X	Y	Bias	Specimen	X	Y	Bias	Specimen	X	Y	Bias
S00001	8	8	0	S00015	8	8	0	S00029	5	5	0
S00002	3	3	0	S00016	132	128	-4	S00030	12	12	0
S00003	11	10	-1	S00017	7	6	-1	S00031	11	11	0
S00004	12	12	0	S00018	5	5	0	S00032	3	3	0
S00005	8	8	0	S00019	6	6	0	S00033	9	9	0
S00006	4	4	0	S00020	6	5	-1	S00034	5	5	0
S00007	6	6	0	S00021	4	4	0	S00035	4	5	1
S00008	12	12	0	S00022	9	9	0	S00036	6	6	0
S00009	14	9	-5	S00023	3	4	1	S00037	11	10	-1
S00010	5	5	0	S00024	12	12	0	S00038	8	8	0
S00011	69	67	-2	S00025	6	6	0	S00039	8	8	0
S00012	4	4	0	S00026	7	7	0	S00040	6	6	0
S00013	17	17	0	S00027	5	5	0				
S00014	191	180	-11	S00028	15	15	0				

Values with an "X" were excluded from the calculations.