# CRYO.S BIOBANKING TUBES

## **Gebrauchsanweisung Instructions For Use Mode d'emploi Istruzioni per l'uso Használati útmutató**

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E-MAIL [info@de.gbo.com](mailto:info@de.gbo.com)

Greiner Bio-One is a global player.

Find the contact details of your local partner on our website.

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\* X in Order no. for tubes with coloured cap: x = 6 - green; x = 5 - yellow; x = 4 - blue; x = 3 - red; x = 8 - pink; x = 9 - brown; x = 7 - black.

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**greiner**  
BIO-ONE

CRYO.S BIOBANKING TUBES, 600 µL, 2D CODES, BULK, WITH SCREW CAP, STERILE, 192 PCS./BAG, 1 MANUAL CAPPING TOOL

**LOT E1234567**  
2026-04-07

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**STERILE** **CE** **UKCA** **IVD**

**HIBC** **GSI UDI**  
(01) 3403964300515  
(10) E1234567  
(17) 260407

**FOR LABORATORY USE ONLY**

**CH REP** **UK REP**

Greiner Bio-One GmbH  
Maybachstr. 2  
72636 Frickenhausen  
Germany

Greiner Bio-One  
Vaudeville Schweiz GmbH  
St. Leonhardstrasse 39  
8000 St. Gallen  
Switzerland

Greiner Bio-One Ltd.  
Brunel Way  
Stourminster Business Park  
GL 10 3SX, Stonehouse  
United Kingdom

Made in Germany Rev.03 - 1/02/23

**REF 132263**  
CRYO.S BIOBANKING TUBES, 600 µL, 2D CODES, BULK, WITH SCREW CAP, STERILE, 192 PCS./BAG, 1 MANUAL CAPPING TOOL

**greiner**  
BIO-ONE

CRYO.S BIOBANKING TUBES, 600 µL, 2D CODES, BULK, WITH SCREW CAP, STERILE, 192 PCS./BAG, 1 MANUAL CAPPING TOOL

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- (GB) **Manual tube opening:** If tubes with pre-installed screw caps are used, screw caps must be removed before pipetting of sample into tubes. Alternatively un-capped tubes (Order nos. 976570, 977570, 978570) may be used in combination with separately available screw caps (Order nos. 385270, 38527X). Remove the rack lid by firmly pressing thumb and forefinger against opposing latches and moving the rack lid upwards.

Fully insert the hexagonal portion of the blue capping/decapping tool (contained in each box) into the central recess of the screw cap. Turn the decapping tool counter-clockwise until the screw cap gets loose. The tube will interlock with rotation stoppers of the rack during this step. The decapping tool with the mounted screw cap may be laid down on an even, clean surface until re-closure of the tube. Alternatively, screw caps may be stored in a clean 96-well U-bottom plate (Order no. 650161). After all pipetting steps from or into the tubes have been completed, screw caps may be re-applied using the blue manual capping/decapping tool and stripped off the tool by pressing the forefinger against the upper surface of the screw cap.

**Machine-assisted Opening and Closing of Tubes:** Cryo.s Biobanking tubes are designed for the opening and closing with semi-automated and automated devices such as the Cryo.s 8-Channel Handheld Decapper from Greiner Bio-One (Order no. 852070) and the Labelite DeCapper from Hamilton Storage Technologies, respectively. 96-well U-bottom plates (Order no. 650161) may be used as cap carriers for the intermediate storage of screw caps in combination with these devices. Torque values of 5-7 Ncm were found suitable for tube closure.

**Scanning of Racks and Tubes:** Cryo.s Biobanking tubes carry a datamatrix code surrounded by a corresponding human readable text on their bottom for tube identification. The corresponding racks carry a linear barcode on their side and a datamatrix code on their bottom; both these codes encode the same rack identifier and may be used for rack identification. For rack scanning use the Cryo.s Rack Scanner (Order no. 849050 [US edition] or 849070 [with international AC adapter plug set]) or other rack scanners capable of scanning racks following the SLAS plate format standard.

**Sample Shipment:** Cryo.s Biobanking tubes may be used as a primary receptacle for the shipping of biological material, including air shipment. Potentially hazardous biological materials must be packaged to withstand content leakage, shocks, temperature changes, pressure changes, and other conditions that can occur during transport. It is important to follow the latest recommendations of the International Air Transporters Association (IATA) including triple packaging and correct labeling. Cryo.s Biobanking tubes alone are not sufficient as sample containers for the shipment of hazardous material, but must be used in combination with suited secondary and tertiary packaging materials.

**Centrifugation:** Cryo.s Biobanking tubes may be centrifuged in the corresponding racks in suitable, precisely fitting swinging bucket rotors with anti-slide pat. Do not centrifuge Cryo.s Biobanking tubes in Racks at more than 3.000 g.

**Maximum working volumes:** Maximum working volumes must not be exceeded. They are: 235  $\mu$ L liquid for all 300  $\mu$ L tube versions, 580  $\mu$ L liquid for all 600  $\mu$ L tube versions, 975  $\mu$ L liquid for all 1000  $\mu$ L tube versions.

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GB

During all steps of work with Cryo.s Biobanking tubes wear laboratory coat and gloves, thus fully covering your arm and hand.

Wear goggles during all steps of work with Cryo.s Biobanking tubes.

Do not use tubes and racks whose primary packaging is damaged, as this can lead to contamination. Please inspect the primary packaging visually before use!

Do not use tubes whose screw cap got loose!

When working with potentially hazardous biological materials apply a laboratory practice with maximum safety precautions that is customary according to the state of the art.

Never submerge Cryo.s™ Biobanking tubes into liquid nitrogen. Cryo.s Biobanking tubes are suited for storage in the gas phase above liquid nitrogen or in freezers only.

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UK REP

Greiner Bio-One Ltd.  
Brunel Way, Stroudwater Business Park  
GL10 3SX, U.K. Stonehouse, Gloucestershire  
United Kingdom

making a difference

[www.gbo.com](http://www.gbo.com)

**GREINER BIO-ONE GMBH**  
FRICKENHAUSEN, GERMANY

**PHONE** +49 7022 948-0  
**FAX** +49 7022 948-514  
**E-MAIL** [info@de.gbo.com](mailto:info@de.gbo.com)

**GREINER BIO-ONE IS A GLOBAL PLAYER.**  
FIND THE CONTACT DETAILS OF YOUR  
LOCAL PARTNER ON OUR WEBSITE.

Devices of Greiner Bio-One are to be used by properly qualified persons only in accordance with the relevant Instructions for Use (IFU), where applicable.

For more information contact your local Greiner Bio-One sales representative or visit our website ([www.gbo.com](http://www.gbo.com)).

All information is provided without guarantee despite careful processing.

Any liability, warranty or guarantee of Greiner Bio-One GmbH is excluded.

All rights, errors and changes are reserved. If not stated otherwise, Greiner Bio-One GmbH has all copyrights and/or other (user-)rights in this documents, in particular to signs such as the mentioned (word-picture-)brands and logos. Any use, duplication or any other use of the rights of Greiner Bio-One GmbH is expressly prohibited.

**Media owner:** Greiner Bio-One GmbH / Represented by Managing Director Bernd Klingel. The company is registered in the Commercial Register at the first instance court in Stuttgart, HRB 224604 / VAT Number: DE812585719.

  
**greiner**  
BIO-ONE

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Greiner Bio-One Vacuette Schweiz GmbH