# CLiP® Automatic safety I.V. catheter

## Instructions for use:

EN	CLiP® Automatic safety I.V. catheter	PT	CLiP® Cateter I.V. de segurança automática
DA	CLiP® automatisk sikkerheds PVK	EL	Αυτόματος ενδοφλέβιος (I.V.) καθετήρας ασφαλείας CLiP®
NO	CLiP® automatisk I.V. sikkerhetskanyle	SR	CLiP® ИВ канила са аутоматском заштитом
SV	CLiP® I.V. kateter med automatisk säkerhetsmekanism	HR	CLiP® Automatic (automatski) sigurnosni intravenozni kateter
FI	Automaattinen laskimon CLiP®-turvakatetri	LT	CLiP® IV kateteris su automatine apsauga
DE	CLiP® Automatischer intravenöser Sicherheitskatheter	TR	CLiP® Otomatik emniyetli I.V. kateteri
NL	CLiP® Automatische I.V. veiligheidskatheter	RO	Cateter I.V. de siguranță automat CLiP®
FR	Cathéter I.V. à mise en sécurité automatique CLiP®	PL	Cewnik dożylny z automatycznym zabezpieczeniem CLiP®
IT	CLiP® Catetere venoso di sicurezza automatico	CS	CLiP® automatický bezpečnostní intravenózní (nitrožilní) katétr
ES	Catéter I.V. de seguridad automático CLiP®		

CLiP® Winged, CLiP® Ported and CLiP® Neo are automatic safety intravenous (I.V.) catheters. The products have a built in safety mechanism that encapsulates the tip of the used needle when extracted from the catheter. This prevents accidental needle stick injury provided that general effective and safe working procedures and precautions are maintained during its use and disposal. The products feature stabilization wings for secure fixation.

#### Intended use

Intravenous/Intravascular access for short term peripheral cannulation (up to 30 days). The use of this product is restricted to qualified healthcare professionals.

#### Indications

- Infusion of I.V. solutions including blood and fluid of similar viscosity.
- Intermittent intravenous drug administration.
- Sample blood
- The products withstand use with power injectors rated for a maximum pressure of 21 bar (305 psi).

#### Contraindications

Product should not be used in patients with known hypersensitivity to any of the materials used.

#### Materials used

Polypropylene (PP), Polyethylene (PE), Polycarbonate (PC), Polyoxymethylene (POM), silicone, stainless steel. Catheter with radiopaque lines: Fluorinated Ethylene Propylene (FEP) or Polyurethane (PUR), see package labelling. This product is not made with natural rubber latex. This product is PVC and DEHP free. This product is sterilized with Ethylene Oxide gas.

#### Instructions for use

- Carefully select and aseptically prepare the site.
- Select suitable I.V. catheter size. Make sure that packaging is not damaged and check the expiry date.
- Remove I.V. catheter from sterile packaging.
- Grip the I.V. catheter, remove the needle cover and make sure that the needle and the catheter tip are properly aligned.
- Perform venipuncture and check for blood indication in flashback chamber (22-14 G) or along the catheter length (26-24 G).
- Advance the catheter, ensure the tip is in the vein. Advance further, hold the catheter hub for stabilization and simultaneously withdraw the needle completely. Always pull the needle straight back while pressing the vein just above the tip of the catheter to prevent spillage of blood. Confirm that the safety mechanism encapsulate the needle tip.
- Discard the needle in approved sharps container.
- Connect I.V. infusion line, connector or cap with a 6% (Luer) taper connection.
- CLiP® Winged and CLiP® Neo: Remove the wing holder.
- Secure the I.V. catheter to ensure clear fluid path and cover the insertion site with a sterile dressing.
- CLiP® Ported: The port can be used for injections using a needleless syringe. Close the port cap after use.
- Perform routine monitoring, venipuncture site care and change intervals according to local procedures. Advise patient to inform a health care professional of any pain, leakage or discomfort from the insertion site.

#### Caution

- The product shall be used immediately after opening the packaging to ensure product sterility.
- I.V. catheters placed near flexing joint and/or poorly fixated can cause kinking of catheter.
- Do not bend the needle before use, during advancement or removal. Bent needle can damage the catheter and/or the safety mechanism functionality.
- Do not re-insert a partially or completely withdrawn needle as this could damage the catheter.
- Covering (e.g. with finger or dressing) the end of the catheter hub during withdrawal could prevent safety mechanism activation.
- For single use only, re-use can cause cross-contamination.
- Prior use with power injector ensure product functionality.
- Do not use scissors or other sharp objects when removing fixation dressing as it may damage the catheter.

#### Label information

The numerals marked next to symbol for EO sterilization on single unit label are the identifier for cavity in packing die.

The specification numbers printed on packaging and labels are graphic reference numbers.

"HR/DEVICE/MFG/MD/2019/000091" Code Number issued by Indian Competent Authority for export labelling.

For symbol explanations see separate table in this instruction for use.

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