

VACUETTE® Disposable Tourniquet Instructions for use





Intended Use

The **VACUETTE®** Disposable Tourniquet is used in routine blood collection to apply pressure to the arm to make the veins more prominent and easier to feel and find. The product is to be used by appropriately trained healthcare professionals only in accordance with these instructions.

Product Description

The **VACUETTE®** Disposable Tourniquet is made of synthetic rubber. It is not made with natural rubber latex to prevent allergies and skin irritation. The pre-perforated parts can be easily torn off as required.

Handling

			
1 Place the tourniquet 7.5-10 cm above the puncture site.	2 Wrap the VACUETTE® Disposable Tourniquet around the arm, crossing the strands over each other	3 Tuck the top strand under the crossed tourniquet to create a loop, without pulling the end through.	4 When the blood flows into the first tube, loosen the tourniquet by pulling the right end in an upward motion.

If a tourniquet has been applied for longer than 1 minute before accessing the vein, it should be released and reapplied after 2 minutes before venipuncture is performed.¹

Precautions/Cautions

- Do not apply high pressures that may impair arterial blood flow to the extremity.
- Release tourniquet as soon as blood starts to flow into the first tube.
- Too long and too high stasis may contribute to haemolysis.
- Do not use on broken skin.
- Dispose tourniquet immediately after use and do not re-use.
- Reuse of the disposable tourniquet can lead to illness for the patient due to cross-contamination.

Only applicable for member states of the European Union: Should any serious incidents occur in relation to the product, these must be reported to the manufacturer and the competent authority in the member state, in which the user/patient is established.







Storage

Storage temperature: up to 25°C. Avoid exposure to direct sunlight.

Disposal

Dispose of according to your facility's policy.

Label information

	Manufacturer		Keep away from sunlight	MD	Medical device
	Use-by date		Upper limit of temperature	CH REP	Authorized representative in Switzerland
LOT	Batch code		Do not re-use	UK REP	UK responsible person
REF	Catalogue number		Consult instructions for use		

Literature:

¹CLSI. *Collection of Diagnostic Venous Blood Specimens*. 7th ed. CLSI standard GP41. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.



Greiner Bio-One GmbH
Bad Haller Str. 32
4550 Kremsmünster
Austria

Production location:

Manufacturing facilities of Surge Industrial Co., Ltd., 4th Fl. No. 258,
Sec. 3, Beishen Road.,
Shenkeng District, New Taipei City 22204, Taiwan,
Distributed by Greiner Bio-One GmbH, Austria

www.gbo.com/preanalytics
office@at.gbo.com
Phone +43 7583 6791



Greiner Bio-One Vacuette Schweiz GmbH
St. Leonhardstraße 39
9000 St. Gallen, Switzerland



Greiner Bio-One Ltd, Brunel Way, Stroudwater Business
Park, GL10 3SX U.K. Stonehouse, Gloucestershire,
United Kingdom

