

Your **Power** for Health

# THE GENTLE TOUCH

The next evolution of safety.

**VACUETTE® EVOPROTECT**

SAFETY Blood Collection Set



  
**greiner**  
BIO-ONE

# VACUETTE® EVOPROTECT

## SAFETY Blood Collection Set

**EVOPROTECT** is the next evolution in blood collection. The semi-automatic safety mechanism protects the user from the risk of needlestick injuries and makes the daily task of blood collection easier. **VACUETTE® EVOPROTECT** provides for gentle collection and safety from injury because we believe that every needlestick is one too many.

The Safety Blood Collection Set meets current safety recommendations and requirements while providing for safe blood collection and infusion with the same set.

### User-friendly with an optimized design

The **VACUETTE® EVOPROTECT** is ergonomically designed and intuitive to use. It has a winged needle and is designed for one-handed activation. The potential for needlestick injuries is minimized by the activation of the safety mechanism while the needle is still in the vein.



### 1 MAXIMUM SAFETY

The **VACUETTE® EVOPROTECT** Safety Blood Collection Set has a semi-automatic safety mechanism. The triggers are located on the side of the product, facilitating one-handed activation. The needle is then automatically retracted and, at the same time, safely and irreversibly enclosed in the safety shield. A clearly audible click confirms that the safety mechanism is successfully activated.

### 2 DESIGNED WITH PATIENT COMFORT IN MIND

The ultra sharp, triple bevel needle tip ensures optimal patient comfort. A flash of blood can be seen between the needle and the tubing if venipuncture is successful. The tubing is designed to ensure the needle stays as still as possible in the vein. The location of the safety mechanism triggers, on the side of the device, minimizes the risk of the needle moving in the patient's vein during activation.

### 3 VERSATILE USE

The same **VACUETTE® EVOPROTECT** Safety Blood Collection Set can be used for drawing blood and for short-term infusion. Using the same set for blood collection and infusion means a second puncture is not needed.

# THE GENTLE TOUCH IN FOUR STEPS

## WINGED NEEDLE

The material and surface design of the wings provide a simple, safe grip and optimal hold during venipuncture.

## SAFETY MECHANISM


The safety mechanism is activated while still in the patient's vein, offering the user the highest level of safety. The design of the safety mechanism promotes a straight removal of the needle from the vein providing for patient comfort.

## TUBING

The tubing is particularly flexible and not made with natural rubber latex or DEHP.

## NEEDLE

The extra-thin walls of the 21G and 23G needles have a positive effect on flow rate and thus the duration of blood collection/infusion.




## STEP 1

Perform venipuncture (with patient's arm in downward position). Flashback will confirm successful venipuncture, depending on the patient's venous pressure.




## STEP 2

Collect blood according to your facility's procedure. After completion, place gauze over collection site without applying pressure.




## STEP 3

The product is intended for in-vein activation. With one hand, activate the safety mechanism by pressing the blue colored release buttons. Take care not to hinder the automatic retraction.





## STEP 4

The safety mechanism will slide backwards until an audible click is heard. The click is a sign that the safety mechanism has been correctly activated. Dispose of device in a sharps container per facility procedure.



For the complete handling recommendation please refer to the Instructions for Use (item n° 980274) and graphic Instructions for Use (item n° 980275).



## STORAGE

4-25°C  
(avoid direct sunlight)

## SHELF LIFE

36 months  
(from day of production)

## PACKAGING

single-packed and sterile

## EASILY-DISPOSABLE

The **VACUETTE® EVOPROTECT** Safety Blood Collection Set packaging is made of one single material, so there is no need to separate paper and plastic for disposal. The soft blister package is also easy to compress, thus reducing waste.

# ITEM OVERVIEW

ITEM NUMBER	DESCRIPTION	GAUGE	NEEDLE LENGTH	TUBING LENGTH	PIECES PER SHELF BOX	PIECES PER CARTON
VACUETTE® EVOPROTECT SAFETY Blood Collection/Infusion Set						
450120	VACUETTE® EVOPROTECT SAFETY Blood Collection/Infusion Set	21G (0.80mm)	19mm 3/4"	7 ½" 19cm	60	1080
450121	VACUETTE® EVOPROTECT SAFETY Blood Collection/Infusion Set	23G (0.60mm)	19mm 3/4"	7 ½" 19cm	60	1080
VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Luer Adapter						
450127	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Luer Adapter	21G (0.80mm)	19mm 3/4"	7 ½" 19cm	60	1080
450128	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Luer Adapter	21G (0.80mm)	19mm 3/4"	12" 30cm	60	1080
450129	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Luer Adapter	23G (0.60mm)	19mm 3/4"	7 ½" 19cm	60	1080
450130	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Luer Adapter	23G (0.60mm)	19mm 3/4"	12" 30cm	60	1080
450132	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Luer Adapter	25G (0.50mm)	19mm 3/4"	12" 30cm	60	1080
VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Holder						
450122	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Holder	21G (0.80mm)	19mm 3/4"	7 ½" 19cm	30	540
450123	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Holder	21G (0.80mm)	19mm 3/4"	12" 30cm	30	540
450125	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Holder	23G (0.60mm)	19mm 3/4"	7 ½" 19cm	30	540
450126	VACUETTE® EVOPROTECT SAFETY Blood Collection Set + Holder	23G (0.60mm)	19mm 3/4"	12" 30cm	30	540

**Austria (Main Office)**  
Greiner Bio-One GmbH  
Phone +43 7583 6791-0  
Fax +43 7583 6318  
E-Mail office@at.gbo.com

**Brazil**  
Greiner Bio-One Brasil  
Phone +55 19 3468-9600  
Fax +55 19 3468-9621  
E-Mail office@br.gbo.com

**China**  
Greiner Bio-One Suns Co., Ltd.  
Phone +86 10 83 55 19 91  
Fax +86 10 63 56 69 00  
E-Mail office@cn.gbo.com

**France**  
Greiner Bio-One SAS  
Phone +33 1 69 86 25 25  
Fax +33 1 69 86 25 35  
E-Mail accueil@gbo.com

**Germany**  
Greiner Bio-One GmbH/  
Preanalytics  
Phone +49 7022 948-0  
Fax +49 7022 948-514  
E-Mail office@de.gbo.com

**Hungary**  
Greiner Bio-One Hungary Kft.  
Phone +36 96 213 088  
Fax +36 96 213 198  
E-Mail office@hu.gbo.com

**India**  
Greiner Bio-One  
INDIA Pvt., Ltd.  
Phone +91 120 375 9291  
E-Mail info@gboindia.com

**Italy**  
Greiner Bio-One Italia S.r.l.  
Phone +39 02 9438 3340  
Fax +39 02 9438 3366  
E-Mail office@it.gbo.com

**Japan**  
Greiner Bio-One Co. Ltd.  
Phone +81 3 3505 8875  
Fax +81 3 3505 8945  
E-Mail info@jp.gbo.com

**Netherlands**  
Greiner Bio-One B.V.  
Phone +31 172 4209 00  
Fax +31 172 4438 01  
E-Mail info@nl.gbo.com

**Portugal**  
VACUETTE Portugal S.A.  
Phone +351 252 647 721  
Fax +351 252 647 722  
E-Mail info@vacuette.pt

**Spain**  
Greiner Bio One España, S.A.U.  
Phone +34 91 652 77 07  
Fax +34 91 652 33 35  
E-Mail info@es.gbo.com

**Switzerland**  
Greiner Bio-One VACUETTE  
Schweiz GmbH  
Phone +41 7 12 28 55 22  
Fax +41 7 12 28 55 21  
E-Mail office@ch.gbo.com

**Thailand**  
Greiner Bio-One Thailand Ltd.  
Phone +66 38 4656 33  
Fax +66 38 4656 36  
E-Mail office@th.gbo.com

**Turkey**  
Greiner Bio-One Turkey  
Phone +90 216 576 6004  
Fax +90 216 576 6005  
E-Mail office@tr.gbo.com

**United Kingdom**  
Greiner Bio-One Ltd.  
Phone +44 1453 8252 55  
Fax +44 1453 8262 66  
E-Mail info@uk.gbo.com

**USA**  
Greiner Bio-One  
North America Inc.  
Phone +1 704 261-7800  
Fax +1 704 261-7899  
E-Mail office@us.gbo.com