



No conversion factor


Unlike additives added in liquid form, the powder additive in the **VACUETTE®** FC Mix tube has no dilution effect. There is therefore no need to take a conversion factor into consideration.

Inverting ten times ensures that the tube additive is completely dissolved and therefore well mixed with the sample.

Tried-and-tested vacuum technology

Greiner Bio-One's tried-and-tested vacuum technology is used for the new **VACUETTE®** FC Mix tube. The shatter-proof tube is made of polyethylene terephthalate (PET). PET is important for the stability of the vacuum. The safety cap is particularly easy to open and allows for hygienic working. The transparent plastic label provides an optimum view of the tube contents.

Greiner Bio-One provides the following centrifugation recommendation for collecting plasma in FC Mix Tubes:



Your **Power** for Health



Austria (Main Office)
Greiner Bio-One GmbH
Phone +43 7583 6791-0
Fax +43 7583 6318
E-Mail office@at.gbo.com

Brazil
Greiner Bio-One Brasil
Phone +55 19 3468-9600
Fax +55 19 3468-9621
E-Mail office@br.gbo.com

China
Greiner Bio-One Suns Co., Ltd.
Phone +86 10 83 55 19 91
Fax +86 10 63 56 69 00
E-Mail office@cn.gbo.com

France
Greiner Bio-One SAS
Phone +33 1 69 86 25 25
Fax +33 1 69 86 25 35
E-Mail office@fr.gbo.com

Germany
Greiner Bio-One GmbH/Preanalytics
Phone +49 7022 948-0
Fax +49 7022 948-514
E-Mail office@de.gbo.com

For further information, please visit our website www.gbo.com/preanalytics or contact us.

Hungary
Greiner Bio-One Hungary Kft.
Phone +36 96 213 088
Fax +36 96 213 198
E-Mail office@hu.gbo.com

India
Greiner Bio-One INDIA Pvt., Ltd.
Phone +91 120 375 9291
E-Mail info@gboindia.com

Italy
Greiner Bio-One Italia S.r.l.
Phone +39 02 9438 3340
Fax +39 02 9438 3366
E-Mail office@it.gbo.com

Japan
Greiner Bio-One Co. Ltd.
Phone +81 3 3505 8875
Fax +81 3 3505 8945
E-Mail info@jp.gbo.com

Netherlands
Greiner Bio-One B.V.
Phone +31 172 4209 00
Fax +31 172 4438 01
E-Mail info@nl.gbo.com

Portugal
VACUETTE Portugal S.A.
Phone +351 252 647 721
Fax +351 252 647 722
E-Mail info@vacuette.pt

Spain
Greiner Bio One España, S.A.U.
Phone +34 91 652 77 07
Fax +34 91 652 33 35
E-Mail info@es.gbo.com

Switzerland
Greiner Bio-One VACUETTE Schweiz GmbH
Phone +41 7 12 28 55 22
Fax +41 7 12 28 55 21
E-Mail office@ch.gbo.com

Thailand
Greiner Bio-One Thailand Ltd.
Phone +66 38 4656 33
Fax +66 38 4656 36
E-Mail office@th.gbo.com

Turkey
Greiner Bio-One Turkey
Phone +90 216 576 6004
Fax +90 216 576 6005
E-Mail office@tr.gbo.com

United Kingdom
Greiner Bio-One Ltd.
Phone +44 1453 8252 55
Fax +44 1453 8262 66
E-Mail info@uk.gbo.com

USA
Greiner Bio-One North America Inc.
Phone +1 704 261-7800
Fax +1 704 261-7899
E-Mail info@us.gbo.com



980295 rev.01 - 02.2018 e

Your **Power** for Health



VACUETTE® FC Mix Tube

Glucose stabilisation
right from the beginning

VACUETTE® FC Mix tube

Glucose stabilisation right from the beginning

The breakdown of glucose (glycolysis) in venous blood samples is of great significance in pre-analytics, particularly in relation to the diagnosis of diabetes mellitus and gestational diabetes. Various guidelines (American Association of Clinical Chemistry, American Diabetes Association, German Diabetes Society [DDG] and German Society for Gynaecology and Obstetrics [DGGG]) deal with this matter in detail^{1,2}.

Greiner Bio-One now has a solution in the form of the **VACUETTE® FC Mix tube**.

According to the guidelines from the DDG*/DGGG**, the citrate fluoride additive in the tube stabilises the *in vivo* concentration of glucose in the sample¹.

The advantages are clear:

- ✓ **Stabilisation immediately after collection for 48 hours**
- ✓ **Based on the *in vivo* value (almost 100%)**
- ✓ **Avoids false negative diagnoses of diabetes patients**
- ✓ **Stabilisation allows for longer transport and storage times**

Effective glycolysis inhibition for precise determination of the *in vivo* blood sugar content



The unique additive mixture is what makes the difference:

Citrate/citric acid buffer ensures quick stabilisation

Glycolysis depends on the pH value. It is catalysed by the enzymes hexokinase and phosphofructokinase. When stored between + 4°C and room temperature, the enzymes, and consequently also glycolysis, are suppressed and the blood sugar is therefore constantly held at the *in vivo* value^{3,4}.

Inhibition via sodium fluoride

In order to extend the inhibition to 48 hours, the tube contains a sodium fluoride additive.

After correct inversion FC Mix tubes can be stored for up to 24 hours at room temperature without centrifugation

Should the tubes be expected to be stored longer than 24 hours at room temperature, samples should be immediately centrifuged after blood collection in order to be stored for up to 48 hours at room temperature. Centrifuged aliquots from FC Mix tubes can be stored for up to 48 hours at room temperature. Tubes should be centrifuged as soon as possible.


Cooling of the samples (4-8°C, 39-46°F) is also suitable for 48 hours glucose stabilization.




EDTA as anticoagulant

Na₂-EDTA is included in the mixture as an anticoagulant, which is effective due to complexation of Ca⁺⁺.

Initial glucose level⁵:



Schematic illustration of glycolysis⁶:



Extensive study material on the **VACUETTE® FC Mix tube** can be found in the Download Center on our website www.gbo.com.

References:

- Sirid Griebenow (2016): Suitability of **VACUETTE® FC Mix Blood Collection Tube** for gestational diabetes testing, Whitepaper Hospital Isala, Zwolle (NL)
- Kellerer et al, Gestationsdiabetes mellitus – Praxisleitlinie der DDG und der DGGG, (Gestational diabetes – practice guidelines of the DDG and DGGG) (2011)
- Sacks et al, Guidelines and Recommendations of Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus, Clinical Chemistry 57:6 (2011)
- Yagmur and Van Helden et al, Effektive Glykolyse-Inhibierung im Citratgepufferten venösen Vollblut und Plasma [Effective glycolysis inhibition in citrate-buffered venous full blood and plasma], Lab Med (2012), 36(3):169-177
- Gambino et al, Acidification of Blood Is Superior to Sodium Fluoride Alone as an Inhibitor of Glycolysis, Clinical Chemistry 55:5 (2009)
- Diagnostic Samples: From The patient to the Laboratory, 4th Edition
- Biochemie, Jeremy M. Berg, John L. Tymoczko, Lubert Stryer, 2007

* Deutsche Diabetes Gesellschaft (German Diabetes Society)
 ** Deutsche Gesellschaft für Gynäkologie und Geburtshilfe (German Society for Gynaecology and Obstetrics)

Item Overview

Item No.	Norm. Vol.	Cap Colour	Ring Colour	Tube Size	Description	Packaging Inner	Packaging Outer
FC Mix							
454510	2ml	pink	white	13 x 75	transparent label, non-ridged	50 pcs.	1200 pcs.
454511	2ml	grey	white	13 x 75	transparent label, non-ridged	50 pcs.	1200 pcs.
454513	3ml	pink	black	13 x 75	transparent label, non-ridged	50 pcs.	1200 pcs.
454514	3ml	grey	black	13 x 75	transparent label, non-ridged	50 pcs.	1200 pcs.