

## Puncture Site

**1** **Prewarming** the puncture site before the blood collection for 3 – 5 minutes (with max. 42 °C) can significantly increase the blood flow.

Finger pad:  
**MiniCollect®**  
Safety Lancets (a) or



**MiniCollect®**  
Lancelino Safety Lancets (b)




Heel: Only with  
**MiniCollect®** Safety Lancets




## MiniCollect®

**2** The volume of blood depends on the length of the lancet and is marked with blood drops.


### Safety Lancets with needle

Safety Lancet, lavender, 28 G, penetration depth 1,25 mm 

Safety Lancet, orange, 23 G, penetration depth 2,25 mm 


### Safety Lancets with blade


Safety Lancet, pink, penetration depth 1 mm 

Safety Lancet, green, penetration depth 1,5 mm 

Safety Lancet, blue, penetration depth 2 mm 

### Lancelino Safety Lancets with needle (sterile)


Safety Lancet light green, 30G, penetration depth 1.2 mm 

Safety Lancet lavender, 25G, penetration depth 1.5 mm 

Safety Lancet blue, 21G, penetration depth 1.8 mm 


Safety Lancet green, 21G, penetration depth 2.4 mm 

### Lancelino Safety Lancets with blade (sterile)


Safety Lancet yellow, 0,8 mm, penetration depth 2 mm 

## Handling


**3** Remove the protective cover of the blade or needle



Position the Safety Lancet (a) or Lancelino Safety Lancet (b)



Puncture the skin by pressing button on Safety Lancet (a) or pressing the Lancelino Safety Lancet against the puncture site until it activates (b)




Discard first drop of blood



## Order of draw

**4** We recommend the following Order of draw:



### Note

Always follow your facility's protocol for order of draw.

### Legal Disclaimer

The information contained in this document is the view of the authors and does not imply a formal endorsement on the part of Greiner Bio-One (or the outside organizations that employ the authors). It is provided for information purposes only and does not constitute medical advice nor is it intended or recommended as a substitute for professional medical advice.

The information contained in this document is not a substitute for informed judgment and the reader shall be solely responsible for decisions that may arise in whole or in part from the use of or reliance upon information contained in the document. To the extent permitted by applicable law, neither Greiner Bio-One nor the authors shall be liable for any decision made, action taken or inaction in reliance on the information, statements, evaluations or recommendations set forth in this document.

# MiniCollect®

## Safety Lancets / Lancelino Safety Lancets

### Choose what suits you best