

# CELLdisc™ Accessories



## CELLlevator™ - Easy and secured CELLdisc™ stacking

The CELLlevator™ is a CELLdisc™ accessory facilitating easy and secured stacking of two CELLdisc™ devices to enable best possible utilization of a given space (e.g. in an incubator).

The CELLlevator™ consists of three identical blue polypropylene segments which assembled form a uniform triangle. Centering ribs at the side of the segments recession guide the insertion assuring a simple and faultless assembly process. Thereafter the stacking device is positioned within the circumferential ring of the top plate of the CELLdisc™. Another CELLdisc™ can now be positioned on top of the CELLlevator™.

Ventilation of the lower CELLdisc™ is assured by the notches on the segments whereas planarity of the upper CELLdisc™ is assured by the elevations on the segments.

The maximum loading capacity of the assembled CELLlevator™ is 8 kg which corresponds a filled 40-layer CELLdisc™. Based on the application the CELLlevator™ is delivered unsterile but can be autoclaved if necessary.

<p><b>Key Facts</b></p> <ul style="list-style-type: none"> <li>➤ Easy and secured CELLdisc™ stacking</li> <li>➤ Maximum loading capacity 8 kg</li> <li>➤ Space saving storage</li> <li>➤ Autoclavable (120 °C, 2 bar)*</li> </ul> <p><small>*Autoclaving the CELLlevator™ more than 3 times is not recommended. After each autoclaving cycle the integrity of the segments has to be analyzed by the user.</small></p>		
--	---	--

Feature	Properties
Material	Polypropylene
Max. Loading Capacity	8 kg (assembled CELLlevator™)
Sterility	No
Autoclavability	Yes (120 °C, 2 bar), max. 3 times*
TP-range	4-37 °C

Ordering Information			
Order no.	Product Description	Qty. per bag	Qty. per case
878071	CELLlevator™ - Stacking device for CELLdisc™	1	9

For further information and/or sample ordering please visit our website [www.gbo.com](http://www.gbo.com) or contact us.

Germany (Main office): [info@de.gbo.com](mailto:info@de.gbo.com) | Austria: [office@at.gbo.com](mailto:office@at.gbo.com) | Belgium: [info.be@gbo.com](mailto:info.be@gbo.com) | Brazil: [office@br.gbo.com](mailto:office@br.gbo.com) | China: [info@cn.gbo.com](mailto:info@cn.gbo.com)  
 France: [accueil.France@gbo.com](mailto:accueil.France@gbo.com) | Hungary: [office@hu.gbo.com](mailto:office@hu.gbo.com) | Italy: [office@it.gbo.com](mailto:office@it.gbo.com) | Japan: [info.JP@gbo.com](mailto:info.JP@gbo.com) | Netherlands: [info.nl@gbo.com](mailto:info.nl@gbo.com)  
 Portugal: [info@vacuette.pt](mailto:info@vacuette.pt) | Spain: [info@es.gbo.com](mailto:info@es.gbo.com) | UK: [info.uk@gbo.com](mailto:info.uk@gbo.com) | USA: [office@us.gbo.com](mailto:office@us.gbo.com)

## CELLdisc™ Accessories

### CELLring™ - Levelling ring compensating surface irregularities

CELLring™ assures an exact horizontal position of CELLdisc™ compensating surface irregularities, e.g. of a working bench or an incubator. Thus, cells and media distribute evenly in every single layer of your CELLdisc™.

#### Key Facts

- Ensures exact planar positioning of CELLdisc™
- Guarantees consistent distribution of cells and media in every single layer



#### Ordering Information

Order no.	Product Description	Qty. per bag	Qty. per case
878075	CELLring™ - Levelling ring compensating surface irregularities, e.g. of a working bench or an incubator	1	3

### CELLhandle™ - Gripper for secure & easy transportation

CELLhandle™ facilitates secure and easy transportation of particularly large-sized CELLdisc™ formats as well as convenient emptying of these.

#### Key Facts

- Gripping device for easy lifting and emptying of large-sized CELLdisc™
- Enables single-hand usage



#### Ordering Information

Order no.	Product Description	Qty. per bag	Qty. per case
878074	CELLhandle™ - Gripper for secure and easy transportation	1	1

For further information and/or sample ordering please visit our website [www.gbo.com](http://www.gbo.com) or contact us.

Germany (Main office): [info@de.gbo.com](mailto:info@de.gbo.com) | Austria: [office@at.gbo.com](mailto:office@at.gbo.com) | Belgium: [info.be@gbo.com](mailto:info.be@gbo.com) | Brazil: [office@br.gbo.com](mailto:office@br.gbo.com) | China: [info@cn.gbo.com](mailto:info@cn.gbo.com)  
 France: [accueil.France@gbo.com](mailto:accueil.France@gbo.com) | Hungary: [office@hu.gbo.com](mailto:office@hu.gbo.com) | Italy: [office@it.gbo.com](mailto:office@it.gbo.com) | Japan: [info.JP@gbo.com](mailto:info.JP@gbo.com) | Netherlands: [info.nl@gbo.com](mailto:info.nl@gbo.com)  
 Portugal: [info@vacuette.pt](mailto:info@vacuette.pt) | Spain: [info@es.gbo.com](mailto:info@es.gbo.com) | UK: [info.uk@gbo.com](mailto:info.uk@gbo.com) | USA: [office@us.gbo.com](mailto:office@us.gbo.com)