G

TRANSTUBE®

Transtube® is a unique transport system for the recovery of aerobes and anaerobes incorporating Amies' (modified Stuarts's) liquid medium. **Transtube®** offers significant advantages over contemporary liquid based transport systems particularly as the transport medium is retained in a soft foam pad at the base of each transport tube.

FORMULA:

Potassium chloride	Sodium chloride
Sodium Thioglycollate	Deionised water
Magnesium chloride	
Di sodium hydrogen phosphate	
Potassium d, hydrogen phosphate	
pH 7.3 ± 0.2	

Material safety information

Transtube plastic components do not contain latex or PVC.

INSTRUCTIONS FOR USE

1. Peel back pouch until applicator and Transtube caps are exposed.
2. Remove Transtube cap between forefinger and thumb and discard.
3. Remove applicator and take sample.
4. Place applicator into Transtube and push down firmly until the Bell cap reaches line.
5. Squeeze bottom of tube to release media
6. Fill in patients details.
7. Transport to laboratory immediately.

FINAL PRODUCT QC

A representative sample of Transtubes are taken at random from every batch and used for testing and reference.

Samples are tested to ensure that the minimal inocula of a spectrum of organisms are recovered after a specified duration of time. Suspensions of the following organisms are absorbed on to the swab of each product sample and placed in its respective medium.

Neisseria gonorrhoeae N.C.T.C. 8375
Haemophilus influenzae N.C.T.C. 8143
Bacteroides fragilis N.C.T.C. 9343
Bordetella pertussis N.C.T.C 0739

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