

Lot/Batch Number List GREEN

Item No	Prod	Code	Additive	Volume	Size	Cap Type	Cap Type	Cap color	Ring Color	Austria	USA	Brazil	Thailand
454008	AAAG	LH	Lithium Heparin Sep	4 mL	13 x 75	NR	pull	Green / T	Yellow	A1005---			
454029	AAAF	LH	Lithium Heparin	4 mL	13 x 75	NR	pull	Green / T	Black	A10040J			D100500Z
454030	AAAL	NH	Sodium Heparin	4 mL	13 x 75	NR	pull	Green / T	Black	A1005---			
454056	AAAF	LH	Lithium Heparin	4 mL	13 x 75	NR	pull	Green / T	Black	A1005---			
454057	AAAG	LH	Lithium Heparin Sep	4 mL	13 x 75	NR	pull	Green / T	Black	A1005---			
454237	AAAF	LH	Lithium Heparin	2 mL	13 x 75	NR	pull	Green / T	White	A1005---			
454244	AAAF	LH	Lithium Heparin	3 mL	13 x 75	NR	pull	Green / T	Black	A1005---			
454247	AAAG	LH	Lithium Heparin Sep	3 mL	13 x 75	NR	pull	Green / T	Yellow	A1005---			
454285	AAAG	LH	Lithium Heparin Sep	4 mL	13 x 75	NR	pull	Green / T	Yellow	A1005---			
454289	AAAG	LH	Lithium Heparin Sep	4 mL	13 x 75	NR	pull	Green / T	Yellow	A1005---			
454302	AAAL	NH	Sodium Heparin	2 mL	13 x 75	NR	pull	Green / T	White	#NV			
454312	AAAF	LH	Lithium Heparin	3 mL	13 x 75	NR	pull	Green / T	White	A1005---			
454316	AAAL	NH	Sodium Heparin	4 mL	13 x 75	NR	pull	Green / T	Black	#NV			
454423	AAAF	LH	Lithium Heparin	3 mL	13 x 75	NR	pull	Green / T	Black	A1005---			
454434	AAAG	LH	Lithium Heparin Sep	4 mL	13 x 75	NR	pull	Green / T	Yellow	#NV			
454435	AAAG	LH	Lithium Heparin Sep	3 mL	13 x 75	NR	pull	Green / T	Yellow	#NV			
456012	AAB5	K2E	K ₂ EDTA Sep	5 mL	13 x 100	NR	pull	Green / T	Yellow	A1005---			
456025	AABG	LH	Lithium Heparin Sep	5 mL	13 x 100	NR	pull	Green / T	Yellow	A1005---			
456028	AABL	NH	Sodium Heparin	6 mL	13 x 100	NR	pull	Green / T	Black	A1005---			
456087	AABG	LH	Lithium Heparin Sep	5 mL	13 x 100	NR	pull	Green / T	Yellow	A1005---			
456088	AABF	LH	Lithium Heparin	6 mL	13 x 100	NR	pull	Green / T	Black	A1005---			D1004008
456216	AABG	LH	Lithium Heparin Sep	5 mL	13 x 100	NR	pull	Green / T	Yellow	A100417Q			
456236	AABN	K3E	K ₃ EDTA	6 mL	13 x 100	NR	pull	Green / T	Black	A1005---			
456271	AABF	LH	Lithium Heparin	6 mL	13 x 100	NR	pull	Green / T	Black	A1005---			
456272	AABG	LH	Lithium Heparin Sep	4 mL	13 x 100	NR	pull	Green / T	Yellow	#NV			
484022	AAAN	K3E	K ₃ EDTA	4 mL	13 x 75	NR	pull	Green / T	Black	A1005---			
484029	AAAF	LH	Lithium Heparin	4 mL	13 x 75	NR	pull	Green / T	Black	A1005---			
31-112	AABL	NH	Sodium Heparin	6 mL	13 x 100	NR	pull	Green / T	Black	A1005---			
31-212	AABL	NH	Sodium Heparin	6 mL	13 x 100	NR	pull	Green / T	Black	A1005---			
454008RP	AAAG	LH	Lithium Heparin Sep	4 mL	13 x 75	NR	pull	Green / T	Yellow	A1005---			
454029H	AAAF	LH	Lithium Heparin	4 mL	13 x 75	NR	pull	Green / T	Black	#NV			
454030H	AAAL	NH	Sodium Heparin	4 mL	13 x 75	NR	pull	Green / T	Black	#NV			
454237H	AAAF	LH	Lithium Heparin	2 mL	13 x 75	NR	pull	Green / T	White	#NV			
456087RP	AABG	LH	Lithium Heparin Sep	5 mL	13 x 100	NR	pull	Green / T	Yellow	A1005---			
456087SGH	AABG	LH	Lithium Heparin Sep	5 mL	13 x 100	NR	pull	Green / T	Yellow	A1005---			
SA1147	AAAF	LH	Lithium Heparin	4 mL	13 x 75	NR	pull	Green / T	Black	A1005---			

All batch number starting with A1005... are already modified

Batch number legend

A Austria Austria
 YY Year Year
 MM Months Months
 XYZ continuous alphanumeric number starting 001 - ending ZZZ

e.g.

A100300A

A A = Austria, B = US Production, C = Brazil, D = Thailand
 10 Year 2010
 3 March
 00A alphanumeric number
 following number is 00B, previous number is 00: