

Lot/Batch Number List LAVENDER

Item No.	Code	Additive	Volume	Tube Size	Tube Type	Cap Type	Cap Color	Ring Color	Austria	USA	Brazil	Thailand
454052	K2E	K ₂ EDTA	1 mL	13 x 75	NR	pull	Lavender / T	White	A1004---			
31-115	K2E	K ₂ EDTA	2 mL	13 x 75	NR	pull	Lavender / T	White	A1003057			
31-215	K2E	K ₂ EDTA	2 mL	13 x 75	NR	pull	Lavender / T	White	A100305C			
454024	K2E	K ₂ EDTA	2 mL	13 x 75	PREM	screw	Lavender	White	A100302F			
454024H	K2E	K ₂ EDTA	2 mL	13 x 75	NR	pull	Lavender	White	A1004---			
454428	K2E	K ₂ EDTA	2 mL	13 x 75	NR	pull	Lavender	White	A1004---			
454007	K2E	K ₂ EDTA	2 mL	13 x 75	PREM	screw	Lavender	White	A1004---			
454047	K2E	K ₂ EDTA	2 mL	13 x 75	NR	pull	Lavender / T	White	A1004---			
454409	K2E	K ₂ EDTA	2 mL	13 x 75	NR	pull	Lavender / T	White	A1004---			
454020	K2E	K ₂ EDTA	3 mL	13 x 75	PREM	screw	Lavender	Black	A100302C			
454246	K2E	K ₂ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A100302Q			
454246CA	K2E	K ₂ EDTA	3 mL	13 x 75	NR	pull	Lavender	Black	A1004---			
454456	K2E	K ₂ EDTA	3 mL	13 x 75	NR	pull	Lavender	White	A1004---			
454411	K2E	K ₂ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A100306Q			
454457	K2E	K ₂ EDTA	3 mL	13 x 75	NR	pull	Lavender	White	A1004---			
454009	K2E	K ₂ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A1003024			
454023	K2E	K ₂ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A100302D			
454032	K2E	K ₂ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Yellow	A1004---			
454209	K2E	K ₂ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A1002001			
SA1190	K2E	K ₂ EDTA	3,5 mL	13 x 75	NR SW	pull	Lavender / T	Black	A1004---			
454440	K2E	K ₂ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---			
454410	K2E	K ₂ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A10021FL			
454210	K2E	K ₂ EDTA	4 mL	13 x 75	PREM	screw	Lavender	White	A100302J			
454286	K2E	K ₂ EDTA	4 mL	13 x 75	NR	pull	Lavender	White	A1004---			
456023	K2E	K ₂ EDTA	6 mL	13 x 100	PREM	screw	Lavender	Black	A100202X			
456274	K2E	K ₂ EDTA	6 mL	13 x 100	PREM	screw	Lavender	Black	A1004---			
456043	K2E	K ₂ EDTA	6 mL	13 x 100	NR	pull	Lavender / T	Black	A1004---			
455045	K2E	K ₂ EDTA	9 mL	16 x 100	NR	pull	Lavender	Black	A1004---			
456243	K2E	K ₂ EDTA	6 mL	13 x 100	NR	pull	Lavender / T	Black	A1004---			
464034	K3E	K ₃ EDTA	1 mL	13 x 75	PREM	screw	Lavender	White	A1004---			
454034	K3E	K ₃ EDTA	1 mL	13 x 75	PREM	screw	Lavender	White	A100302G			
454403	K3E	K ₃ EDTA	1 mL	13 x 75	PREM	screw	Lavender	White	A1004---			
464087	K3E	K ₃ EDTA	2 mL	13 x 75	PREM	screw	Lavender	White	A1004---			
484500	K3E	K ₃ EDTA	2 mL	13 x 75	NR	pull	Lavender	White	A1004---			
454087	K3E	K ₃ EDTA	2 mL	13 x 75	PREM	screw	Lavender	White	A100302A			
474232	K3E	K ₃ EDTA	3 mL	13 x 75	PREM	screw	Lavender	White	A1004---			
474087	K3E	K ₃ EDTA	2 mL	13 x 75	PREM	screw	Lavender	White	A100302V			
454222	K3E	K ₃ EDTA	2 mL	13 x 75	NR	pull	Lavender / T	White	A1002004			
454458	K3E	K ₃ EDTA	2 mL	13 x 75	NR	pull	Lavender	White	A1004---			
454222H	K3E	K ₃ EDTA	2 mL	13 x 75	NR	pull	Lavender / T	White	A1004---			
454037	K3E	K ₃ EDTA	2 mL	13 x 75	PREM	screw	Lavender	White	A1004---			
454087SGH	K3E	K ₃ EDTA	2 mL	13 x 75	PREM	screw	Lavender	White	A1004---			
454215	K3E	K ₃ EDTA	2 mL	13 x 75	PREM	screw	Lavender	White	A1004---			
454308	K3E	K ₃ EDTA	2 mL	13 x 75	R	screw	Lavender	White	A1004---			
454317	K3E	K ₃ EDTA	2 mL	13 x 75	NR	pull	Lavender / T	White	A1004---			
454015	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Red	A1003020			
474086	K3E	K ₃ EDTA	3 mL	13 x 75	PREM	screw	Lavender	Black	A1003045			
464086	K3E	K ₃ EDTA	3 mL	13 x 75	PREM	screw	Lavender	Black	A10030SO			
484501	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender	Black	A1004---			
454086	K3E	K ₃ EDTA	3 mL	13 x 75	PREM	screw	Lavender	Black	A1003028			
454207	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender	Black	A1004---			
454436	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---			
454217	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1002002			
454459	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender	White	A1004---			
454452	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---			

454217H	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
464217U	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454424	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454417	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Lavender	A1004---
454405	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454232	K3E	K ₃ EDTA	3 mL	13 x 75	PREM	screw	Lavender	White	A1002024
454287	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454460	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender	White	A1004---
454290	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454216	K3E	K ₃ EDTA	3 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
454401	K3E	K ₃ EDTA	3 mL	13 x 75	PREM	screw	Lavender	Red	A1004---
484217	K3E	K ₃ EDTA	3 mL	13 x 75	NR	pull	Lavender / T	Black	A100303G
454103	K3E	K ₃ EDTA	mL (no vac)	13 x 75	PREM	screw	lavender snap cap		A1004---
454003	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Red	A100302T
454425	K3E	K ₃ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454021	K3E	K ₃ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A1002000
454439	K3E	K ₃ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454021H	K3E	K ₃ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454036	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A10021FM
474036	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A10030GH
454448	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
484036U	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
454406	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
454048	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A091103A
454055	K3E	K ₃ EDTA	4 mL	13 x 75	R	screw	Lavender	Black	A1004---
454408	K3E	K ₃ EDTA	4 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
454307	K3E	K ₃ EDTA	4 mL	13 x 75	R	screw	Lavender	Black	A1004---
454036SGH	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
464036	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
464036U	K3E	K ₃ EDTA	4 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
454039	K3E	K ₃ EDTA	4,5 mL	13 x 75	PREM	screw	Lavender	Black	A100302H
454039KWT	K3E	K ₃ EDTA	4,5 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
454223	K3E	K ₃ EDTA	4,5 mL	13 x 75	NR	pull	Lavender / T	Black	A100301Z
454309	K3E	K ₃ EDTA	4,5 mL	13 x 75	R	screw	Lavender	Black	A1004---
456034	K3E	K ₃ EDTA	5 mL	13 x 100	PREM	screw	Lavender	Black	A100304M
456053	K3E	K ₃ EDTA	6 mL	13 x 100	NR	pull	Lavender / T	Black	A1004---
466036	K3E	K ₃ EDTA	6 mL	13 x 100	PREM	screw	Lavender	Black	A1004---
456256	K3E	K ₃ EDTA	6 mL	13 x 100	PREM	screw	Lavender	Black	A1004---
456036	K3E	K ₃ EDTA	6 mL	13 x 100	PREM	screw	Lavender	Black	A100302F
476736	K3E	K ₃ EDTA	6 mL	13 x 100	PREM	screw	Lavender	Black	A1004---
456038	K3E	K ₃ EDTA	6 mL	13 x 100	NR	pull	Lavender / T	Black	A100300R
456038H	K3E	K ₃ EDTA	6 mL	13 x 100	NR	pull	Lavender / T	Black	A1004---
456276	K3E	K ₃ EDTA	6 mL	13 x 100	NR	pull	Lavender / T	Black	A1004---
455226	K3E	K ₃ EDTA	9 mL	16 x 100	NR	pull	Lavender	Black	A1004---
455036	K3E	K ₃ EDTA	9 mL	16 x 100	NR	pull	Lavender	Black	A100300N
455036H	K3E	K ₃ EDTA	9 mL	16 x 100	NR	pull	Lavender	Black	A1004---
465036	K3E	K ₃ EDTA	9 mL	16 x 100	NR	pull	Lavender	Black	A1004---
454035	K3E	K ₃ EDTA Liq	2,5 mL	13 x 75	PREM	screw	Lavender	Black	A1004---
454041	K3E	K ₃ EDTA Liq	3 mL	13 x 75	PREM	screw	Lavender	Black	A10030EA
454042	K3E	K ₃ EDTA Liq	3 mL	13 x 75	NR	pull	Lavender / T	Red	A1004---
454038	K3E	K ₃ EDTA Liq	4 mL	13 x 75	PREM	screw	Lavender	Black	A10030SR
454250	K3E	K ₃ EDTA Liq	4,5 mL	13 x 75	NR	pull	Lavender / T	Black	A1004---
456037	K3E	K ₃ EDTA Liq	6 mL	13 x 100	PREM	screw	Lavender	Black	A10030FE
456050	K3E	K ₃ EDTA Liq	6 mL	13 x 100	NR	pull	Lavender	Yellow	A1004---
456251	K3E	K ₃ EDTA Liq	5 mL	13 x 100	NR	pull	Lavender	Red	A1004---
456209	K3E	K ₃ EDTA Liq	6 mL	13 x 100	NR	pull	Lavender / T	Red	A1004---
454059	K2E	K ₂ EDTA Sep	4 mL	13 x 75	NR	pull	Lavender / T	Yellow	A1004---
454235	K2E	K ₂ EDTA Sep	4 mL	13 x 75	PREM	screw	Lavender	Yellow	A10030ZW
456011	K2E	K ₂ EDTA Sep	5 mL	13 x 100	PREM	screw	Lavender	Yellow	A1004---
456282	K2E	K ₂ EDTA Sep	5 mL	13 x 100	NR	pull	Lavender	Yellow	A1004---
456024	K2E	K ₂ EDTA Sep	4,5 mL	13 x 100	NR	pull	Lavender / T	Yellow	A1004---
455040	K2E	K ₂ EDTA Sep	8 mL	16 x 100	NR	pull	Lavender	Yellow	A1004---
454075	QuantifERON MIT		1 mL	13 x 75	PREM	screw	Lavender	White	A100310N
454414	QuantifERON MIT		1 mL	13 x 75	PREM	screw	Lavender	Yellow	0
370523	-	Snap Cap lavender	-	-	-	-	Lavender	-	

All batch number starting with A1004... are already modified

Batch number legend

A Austria
 YY Year
 MM Months
 XYZ continuous alphanumeric number starting 001 - ending ZZZ

e.g.
 A100300A
 A = Austria, B = US Production, C = Brazil, D = Thailand
 10 Year 2010
 3 March
 00A alphanumeric number
 following number is 00B, previous number is 009