

Lot/Batch Number List LIGHT BLUE

Item No	Prod	Code	Additive	Volume	Size	Cap Type	Cap Type	Cap color	Ring Color	Austria	USA	Brazil	Thailand
454240	AAAA	Z	No Additive	2 mL	13 x 75	NR	pull	Light Blue / T	Black	A100302Y			
454322	AAAR	9NC	Coagulation sodium citrate 3.2%	2 mL	13 x 75	NR SW	pull	Light Blue / T	White	A100300C			
454323	AAAR	9NC	Coagulation sodium citrate 3.2%	2 mL	13 x 75	NR SW	pull	Light Blue / T	White	A1005---			
454324	AAAR	9NC	Coagulation sodium citrate 3.2%	3 mL	13 x 75	NR	pull	Light Blue / T	White	A1003009			
454326	AAAR	9NC	Coagulation sodium citrate 3.2%	4 mL	13 x 75	NR	pull	Light Blue / T	Black	A100300A			
454331	AAAR	9NC	Coagulation sodium citrate 3.2%	3,5 mL	13 x 75	NR SW	pull	Light Blue / T	Black	#NV			
454332	AAAR	9NC	Coagulation sodium citrate 3.2%	3,5 mL	13 x 75	NR SW	pull	Light Blue / T	Black	A100300F			
454333	AAAR	9NC	Coagulation sodium citrate 3.2%	3 mL	13 x 75	NR SW	pull	Light Blue / T	White	A10040N6			
454334	AAAR	9NC	Coagulation sodium citrate 3.2%	3 mL	13 x 75	NR SW	pull	Light Blue / T	Black	A100300I			
454339	AAAR	9NC	Coagulation sodium citrate 3.2%	3 mL	13 x 75	NR SW	pull	Light Blue / T	Black	#NV			
454351	AAAR	9NC	Coagulation sodium citrate 3.2%	3 mL	13 x 75	NR SW	pull	Light Blue / T	Black	#NV			
454352	AAAR	9NC	Coagulation sodium citrate 3.2%	3,5 mL	13 x 75	NR SW	pull	Light Blue / T	White	A1003042			
454360	AAAR	9NC	Coagulation sodium citrate 3.2%	2 mL	13 x 75	NR SW	pull	Light Blue / T	White	A100304C			
454382	AAAS	9NC	Coagulation sodium citrate 3.8%	2 mL	13 x 75	NR SW	pull	Light Blue / T	White	A10031D3			
454386	AAAS	9NC	Coagulation sodium citrate 3.8%	4 mL	13 x 75	NR	pull	Light Blue / T	Black	A10040KJ			
454392	AAAS	9NC	Coagulation sodium citrate 3.8%	3,5 mL	13 x 75	NR SW	pull	Light Blue / T	Black	A100300L			
454393	AAAS	9NC	Coagulation sodium citrate 3.8%	3 mL	13 x 75	NR SW	pull	Light Blue / T	Black	A1005---			
454441	AAAT		CTAD	3 mL	13 x 75	NR SW	pull	Light Blue / T	Yellow	#NV			
484334	AAAR	9NC	Coagulation sodium citrate 3.2%	3 mL	13 x 75	NR SW	pull	Light Blue / T	Black	A100302V			
31-113	AAAS	9NC	Coagulation sodium citrate 3.8%	2 mL	13 x 75	NR SW	pull	Light Blue / T	White	A1003055			
31-114	AAAR	9NC	Coagulation sodium citrate 3.2%	2 mL	13 x 75	NR SW	pull	Light Blue / T	White	A1003056			
31-213	AAAS	9NC	Coagulation sodium citrate 3.8%	2 mL	13 x 75	NR SW	pull	Light Blue / T	White	A100305A			
31-214	AAAR	9NC	Coagulation sodium citrate 3.2%	2 mL	13 x 75	NR SW	pull	Light Blue / T	White	A100305B			
464332U	AAAR	9NC	Coagulation sodium citrate 3.2%	3,5 mL	13 x 75	NR SW	pull	Light Blue / T	Black	A1005---			

484508
SA1190

A1005---
A1005---

All batch number starting with A1005... are already modified

Batch number legend

A Austria Austria
 YY Year Year
 MM Months Months
 XYZ continuous alphanumeric number starting 001 - ending ZZZ

e.g.
 A100300A
 A A = Austria, B = US Production, C = Brazil, D = Thailand
 10 Year 2010
 3 March
 00A alphanumeric number
 following number is 00B, previous number is 009