

Lot/Batch Number List YELLOW

Item No	Prod	Code	Additive	Volume	Size	Cap Type	Cap Type	Cap color	Ring Color	Austria	USA	Brazil	Thailand
454058	AAAD	Z	Serum Sep Clot Activator	4 mL	13 x 75	NR	pull	Yellow / T	Yellow	A1005---			
454205	AAAD	Z	Serum Sep Clot Activator	2,5 mL	13 x 75	NR	pull	Yellow / T	White	A1005---			
454239	AAA7		Urine Culture	4 mL	13 x 75	NR	pull	Yellow / T	Black	A1005---			
454264	AAAD	Z	Serum Sep Clot Activator	4 mL	13 x 75	NR	pull	Yellow / T	Yellow	A1005---			
454280	AAA3		CPDA	4 mL	13 x 75	NR	pull	Yellow / T	Black	A1005---			
456057	AAB3		CPDA	6 mL	13 x 100	NR	pull	Yellow / T	Black	A100417R			
456065	AAB6	Z	Urine No Additive	6,5 mL	13 x 100	NR	pull	Yellow / T	Yellow	A1005---			
456204	AAB2		ACD-B	6 mL	13 x 100	NR	pull	Yellow / T	Black	A1005---			
456238	AABN	K3E	K ₃ EDTA	6 mL	13 x 100	NR	pull	Yellow / T	Black	A1005---			

456200

A1005---

All batch number starting with A1005... are already modified

Batch number legenc

A Austria Austria
 YY Year Year
 MM Months Months
 XYZ continuous alphanumeric number starting 001 - ending ZZZ

e.g.

A100300A

A A = Austria, B = US Production, C = Brazil, D = Thailand
 10 Year 2010
 3 March
 00A alphanumeric number
 following number is 00B, previous number is 009