

Your Power for Health



Your Power for Health




saïnt & baby saïnt

PREMIUM single-use anti-infection tourniquets

Greiner Bio-One Ltd.
Brunel Way
Stroudwater Business Park
Stonehouse
Glos. GL10 3SX
Tel: 01453 825255
Fax:01453 826266
email: sales@uk.gbo.com

VA1203 - printed September 2009



Helping to protect
those most vulnerable
from Healthcare
Associated Infection

VACUETTE®

meet the saint . . .

the new standard in disposable tourniquets

A unique design means that the saint is easy to secure and easy to release, ensuring compliance with local infection control policies. It is ideal for use with the groups highlighted in the **NHS Saving Lives** policy document - Screening for MRSA Colonisation.

High Risk Groups

- Pre-op patients in certain specialities such as Orthopaedics, Cardiothoracic and Neurosurgery
- Emergency Ortho and Trauma admissions especially the elderly and those from nursing or care homes *
- Critical Care (including ICU and HDU)
- Renal Medicine. Patients on renal dialysis have a very high risk of MRSA

Other Specified Groups

- All patients previously known to be MRSA positive
- All elective surgical patients
- Oncology/Chemotherapy in-patients
- Patients admitted from high risk settings especially those admitted from nursing or care homes*

* Some studies have shown 20% of nursing and care home residents are colonised with MRSA

every saint is . . .

Single-use - eliminating risk of patient cross-contamination

Latex-free - providing protection from latex allergy and skin irritation



Sinia

. . and baby saint

for those patients with smaller limbs



Despite its smaller dimensions, baby saint provides all the benefits the saint disposable tourniquet has to offer. Baby saint is ideal for the small limbed patient.

Young or old,
Frail or small,
Baby Saint fits them all !



VACUETTE®

Instructions for Use

The **saint** is only for use by appropriately-trained health professionals, (if in doubt, seek expert advice before use)



1 Simply wrap the **saint** around the arm, passing one end through the slit



2 Hold the slit and pull the notched end firmly through until the desired tension is reached, then release the slit



3 Pull the notched end back on itself and gently release



4 The **saint** will then lock in position



5 When you need to remove it, simply pull the short slitted end and it will release immediately



6 If you need to re apply the **saint**, simply begin again at Step 1

When procedure is complete discard and select a new **saint** for the next patient

Ordering Details

Item No.	Description	Pack Quantity	NPC code
840 054	Saint - Single-use Anti-infection Tourniquet	100	FWJ 018
840 055	Baby Saint - Single-use Anti-infection Tourniquet for patients with smaller limbs	200	FWJ 025