

 GMDN code 10941	TECHNICAL DATA SHEET	309202 en	Date: 02/02/11
		Revision 5	Page: 1/2

Codes: 309202, 409202, 309201, 409201

Screw cap container 25 x 90 mm

2011 catalogue page: 152, 153

Code	Description	Case quantity	Case weight (kg)	Case volume (m ³)	Cases per Pallet
309202	Polystyrene	400	5.30	0.041	48
409202	Polypropylene	400	4.70	0.041	48
309201	Polystyrene, labelled	400	5.50	0.041	48
409201	Polypropylene, labelled	400	4.90	0.041	48

Sanitary Product. IVDD 98 / 79 / CE


SPECIFIC USE, DIMENSIONS AND DESCRIPTION

Codes :	309202 and 309201	409202 and 409201
Material :	CRYSTAL POLYSTYRENE	TRANSPARENT POLYPROPYLENE
Ext. diameter :	25.1 mm.	24.8 mm.
Ext. diameter with cap:	30.8 mm.	30.8 mm.
Total length (C) :	92.5 mm.	91.6 mm.
Maximum volume (D) :	31.7 ml.	31.1 ml.
Recommended volume (E) :	25 ml	25 ml
Graduations :	NO	NO
Withstand centrifugation :	7.000 g	10.000 g

- Conical bottom skirted container suitable for the transport of liquid or solid samples from the collection site up to the laboratory where they will be analysed.
- Each container features a continuous thread and a white ribbed polypropylene cap. Leakproof seal.
- Volume of the conical bottom : 1 ml.

Codes 309201 and 409201: Container especially designed for hospital use.

Labelled container (Minimum order 30.000 units). The label is designed to write the patient's name, the ward and its number, the type of specimen, the time and the date. Dimensions: 50 x 30 mm.



Elaborated	Verified	Approved
