


Technical Product Information



Blood Sampler for pH/blood gas, electrolyte, oximetri and metabolite analysis.

Name of Product:	Blood Sampler for Blood Gas
Volume:	2 ml
Needle:	None
Sterilization:	Gamma Sterilization
Preparation of the Syringe:	Balanced Heparin 50 I.U. per 1 ml blood
Article-No.:	3999599
Shelf Life:	3 years
Packing Unit:	50 pieces per box, 1000 pieces per carton



<u>List of Elements</u>	<u>Material</u>
Body of Device (3)	PP
Luer Adaptor (2)	PP
Closing Cap (1)	PE
Plunger (4)	PP

Manufacturer: SC-Sanguis Counting GmbH
 Im Löchelchen 11
 D-51588 Nümbrecht
 Tel. +49 2293 91290, Fax +49 2293 912950