## **Technical sheet**



## **Lugol, PVP-stabilized Solution**

Code product: 367400

## **Specifications:**

Appearance:

Absorption Maximum in water:

Clear brown solution 354 nm  $\leq \lambda \leq$  357 nm

## **Applications:**

**Bacteriology:** used in Kits Gram-Hücker and Gram-Nicolle stainings for the detection of Gram-positive and Gram-negative bacteria. Lugol, PVP-stabilized Solution forms a complex with Crystal Violet Oxalate that stains Gram-positive bacteria violet.

**Histo-Cytology:** used in May-Grünwald Giemsa staining for panoptic staining of histological sections.

Updated: 25/04/2013