

 DELTALAB	TECHNICAL DATA SHEET	210003 en	Date: 05/03/04
		Revision 5	Page: 1/1

Codes: 210003, 210002, 210004, 210023, 210022, 210024
Fine tip Pasteur pipettes

2003-2004 catalogue page: 92

	Code	Length mm	Capacity ul	Packaging	Case weight	Case volume
900610	210003	104	1,5	8 X400	2.4	0.023
900612	210002	144	5,0	10 X 500	7.4	0.099
900614	210004	153	5,5	8 X 250	2.5	0.040
	210023	104	1,5	500	1.7	0.03
900613	210022	144	5,0	500	1.7	0.03
	210024	153	5,5	500	1.9	0.03

	Code	stem diameter mm	Bulb draw ul	drops ml
900610	210003,210023	3,0	1,0	50
900612 / 900613	210002,210022	5,0	3,3	50
900614	210004,210024	3,0	3,0	50

Sanitary Product. IVDD 98 / 79 / CE

SPECIFIC USE AND DESCRIPTION:

Fine tip pipettes excellent for dispensing micro volumes and small drops (50 drops per ml., i.e. a volume of 0,02 ml./drop).

Material: low density inert polyethylene (LDPE). Low adhesion surface.
 Colour: transparent natural.

Pipettes 210003, 210023, 210004 and 210024 have an **extended fine tip**.
 Pipettes 210002 and 210022 have a **fine tip**.

Can be frozen in liquid nitrogen.

210003/23, 210002/22.- Ideal for Kits-Test, for example HIV, pregnancy; serum separation in Paediatrics mainly, chromatography cup filling, microvials, microplates, electrophoresis gels, E.L.I.S.A., etc...

210004/24.- Ideal for PCR techniques, E.L.I.S.A., DNA samples, protein assays, hybridoma, atomic absorption, and all types of additions in microtechniques, electrophoresis, etc.

210023, 210022, 210024: supplied individually wrapped in a peel-pack sterile by ethylene oxide. On the paper side are printed the product code, batch number and expiry date.

210003, 210002, 210004: supplied in non sterile bags.

A promotional pack featuring our range of pipettes is available, ask our commercial department.

These are not calibrated precision instruments.

Elaborated by Marketing	Authorised by Management