



Specification

GM2 / GM05 1.09204.0000 Giemsa's azur eosin methylene blue solution for microscopy

Specification

Suitability for microscopy (Blood smear)	passes test
Erythrocytes	reddish-brownish
Nuclei	violet
Eosinophilic granules	reddish to red-brown
Neutrophilic granules	red-violet
Lymphocyte cytoplasm	blue
Monocyte plasma	greyblue
Thrombocytes	violet
Basophilic granules	darkviolet
Blood parasite nuclei	lightred

Dr. Karl-August Reiffen

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.