



Product Specification Sheet

VWR TRANSFER PIPETTES, HIGH-PERFORMANCE, GENERAL PURPOSE

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

	Cat. No.	Capacity ml	Drop vol. µl	Suction volume ml	Stem Ø mm	Length mm	Sterile	Packed	PK
PPP1C	612-2846	1,5	40	1	5	87	-	Bulk	1000
	612-2849	1,5	35,7	1	5	123	-	Bulk	1000
	612-2859	23	45,5	10	9	300	-	Bulk	600

Material: LDPE

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.