



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	
GENERAL DESCRIPTION	UTM MLS 3ML PURPLE CAP+FLEXIBLE FLOQSWAB B100
DETAILED DESCRIPTION	1 outer box contain 6 inner boxes each composed of 50 paper-film pouches each containing 1 flexible flocked swab, individually wrapped in paper-film pouch, and 1 tube, with purple screw cap, filled with 3mL of UTM medium.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

Medium pH	7.3 ± 0.2
------------------	-----------

SAMPLING TUBE	CONICAL TUBE	SPECIFICATIONS
		TUBE LENGTH: 101.50 mm ± 0.20 mm TUBE DIAMETER: 16.30 mm ± 0.20 mm INSIDE DIEMETER: 14.10 mm ± 0.20 mm MATERIAL: POLYPROPYLENE

This bulletin (Rev.00-12/07/2023) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

CAP	PURPLE CAP	SPECIFICATIONS
		HEIGHT CAP: 16.30 mm \pm 0.20 mm CAP DIAMETER: 21.50 mm \pm 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE
APPLICATOR	FLEXIBLE FLOCKED SWAB	SPECIFICATIONS
		TOTAL LENGTH: 152.00 mm \pm 2.00 mm STICK DIAMETER: 2.50 mm \pm 0.10 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm FLOCKED DIAMETER: 3.50 mm \pm 0.50 mm BREAK POINT: 98.50 mm \pm 1.00 mm SLIM DIAMETER: 0.90 mm \pm 0.10 mm STICK MATERIAL: POLYSTYRENE FLOCKED MATERIAL: NYLON

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C – DO NOT FREEZE
SHELF LIFE	18 MONTHS

This bulletin (Rev.00-12/07/2023) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.