



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	
GENERAL DESCRIPTION	FECAL 2ML + REGULAR FLOQSWABS BP80 MLS
DETAILED DESCRIPTION	1 outer box contain 10 vipak each composed of 50 film-film pouches each containing 1 flocked regular swab and 1 tube, with green screw cap, filled with 2mL of FECAL medium.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
FECAL	CHLORIDE SALTS SODIUM SALTS PHOSPHATE BUFFER L-CYSTEINE AGAR DISTILLED WATER

Medium pH	7.2 ± 0.3
------------------	-----------

SAMPLING TUBE	TUBE	SPECIFICATIONS
		TOTAL LENGTH: 81.30 mm ± 0.40 mm OUTER TUBE DIAMETER: 13.00 mm ± 0.20 mm INNER TUBE DIAMETER: 11.50 mm ± 0.20 mm MATERIAL: POLYPROPYLENE



TECHNICAL DATA SHEET

CAP	GREEN CAP	SPECIFICATIONS
		CAP HEIGHT: 15.00 mm ± 0.20 mm LOWER CAP DIAMETER: 16.80 mm ± 0.20 mm HIGHER CAP DIAMETER: 16.00 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE
APPLICATOR	FLOCKED REGULAR SWAB	SPECIFICATIONS
		TOTAL LENGTH: 151.00 mm ± 2.00 mm STICK DIAMETER: 2.50 mm ± 0.10 mm BREAKING POINT: 78.50 mm ± 1.00 mm FLOCKED LENGTH: 16.00 mm ± 2.00 mm FLOCKED DIAMETER: 5.50 mm ± 0.50 mm STICK MATERIAL: COPOLYESTER FLOCKED MATERIAL: NYLON

PACKAGING SPECIFICATION

PACKAGING	FILM-FILM POUCH
VI-PAK QUANTITY	50
BOX QTY	500

IMPORTANT INFO

STERILIZATION METHOD	IONIZING IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	15 MONTHS

This bulletin (Rev.01-25/01/2024) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.