

# LIST OF DISINFECTANTS

## VACUETTE® EASYCLEAN

ITEM NO. 840052

### Validated Disinfectants

- Schülke mikrozid® universal wipes premium

### Compatible Disinfectants

- B. Braun Meliseptol® acute Wipes
- B. Braun Meliseptol® sensitive Wipes
- Bode Kohrsolin® FF
- Bode Mikrobac® Tissues
- Clinell® Peracetic acid wipes
- Clinell® universal wipes
- Dr. Schumacher Cleanisept® Wipes Maxi
- Dr. Schumacher Ultrasol OXY® Wipes
- Ecolab Incidin™ Active (powder)
- Ecolab Incidin™ Alcohol Wipe
- Ecolab Incidin™ OxyWipe S
- Ecolab Incidin™ Pro
- Ecolab Incidin™ Rapid
- Ecolab Sekusept™ aktiv (powder)
- Hartmann Bacillol® Sensitive Tissues
- Milton™ Sterilising Tablets
- Sani-Cloth Disinfectant Wipe
- Schülke mikrozid® universal wipes premium
- Schülke perform® (powder)
- SoChlor® DST (1.7g NaDCC)

### Incompatible Disinfectants

- Bode Bacillol® AF

**Note:** Printed copies will not systemically be updated. Valid at the date of issue. Information is subject to change.

## Note on the Disinfectant List

This disinfectant list is intended for use exclusively in conjunction with the methods described in the Instructions for Use (IFU).

The products listed in 'Compatible Disinfectants' have been tested by the manufacturer for material compatibility with the VACUETTE® EasyClean Tourniquet. A complete validation of the reprocessing procedure in combination with the VACUETTE® EasyClean Tourniquet has not been performed. The list serves as an aid to the user in selecting and validating suitable disinfectants.

Responsibility for validating the reprocessing procedure under the respective operating conditions lies with the user or the institution.

The cleaning/disinfecting agents listed are not marketed by Greiner Bio-One; their inclusion is for informational purposes only and does not constitute a recommendation or guarantee.

The use of untested or incompatible disinfectants may lead to premature material fatigue, functional impairment, or product defects.

## Always follow these instructions:

Follow the IFU of the manufacturer of the cleaning agents or disinfectants used (concentration, contact time, temperature, compatibility).

Wear the prescribed personal protective equipment (PPE) during reprocessing of the VACUETTE® EasyClean Tourniquet.

Ensure that no residues of cleaning agents or disinfectants remain on the VACUETTE® EasyClean Tourniquet before use.

Do not use procedures or agents that are not expressly approved in the VACUETTE® EasyClean Tourniquet IFU (e.g., solvents, abrasive agents, non-recommended concentrations).

For technical questions regarding cleaning agents and disinfectants, please contact the respective manufacturer; for questions regarding the VACUETTE® EasyClean Tourniquet, contact us at <https://www.gbo.com/en-int/contact>.

## Disclaimer

To the extent permitted by law, Greiner Bio-One, in its role as distributor of the VACUETTE® EasyClean Tourniquet, assumes no warranty, representation, or liability for damage to materials, equipment, or persons resulting from the use of this list, improper reprocessing, or the use of disinfectants contrary to their IFU.

The cleaning agents and disinfectants listed are third-party products; their respective manufacturers remain responsible for safety, suitability, and compliance with their IFU.

Mandatory statutory rights and liabilities (especially under product liability law) remain unaffected. The applicable references are always the VACUETTE® EasyClean Tourniquet IFU, the IFU of the cleaning agent and disinfectant manufacturers, and relevant legal and normative requirements.