
 greiner bio-one Page 1 / 4	<b>SAFETY DATA SHEET</b> acc. to ISO 11014, 29 CFR 1910.1200	 <b>greiner bio-one</b> A AN SOP 04.03.02-042 Rev.01 Valid from: Sept 12, 2016
	<b>VACUETTE® 9NC Coagulation Sodium Citrate 3,2%          Blood Collection Tube</b>	

## SECTION 1: PRODUCT AND COMPANY IDENTIFICATION

- **Product Name:** VACUETTE® 9NC Coagulation Sodium Citrate 3,2% Blood Collection Tube
- **Manufacturer/Supplier:**

### AUSTRIA

Greiner Bio-One GmbH  
 Bad Haller Strasse 32  
 4550 Kremsmünster  
 Austria  
 Tel: (+43) 7583 6791-0  
 Fax: (+43) 7583 6791-1114  
 Email: [office@at.gbo.com](mailto:office@at.gbo.com)  
 Emergency phone number: (+43) 7583 6791-0

### USA

Greiner Bio-One North America Inc.  
 4238 Capital Drive  
 Monroe, NC 28110  
 USA  
 Tel: (+1) 888-286-3883  
 Fax: (+1) 800-726-0052  
 Email: [info@us.gbo.com](mailto:info@us.gbo.com)  
 Emergency phone number: (+1) 888-286-3883

### BRASIL

Greiner Bio-One Brasil  
 Produtos Médicos Hospitalares Ltda.  
 Av. Affonso Pansan no. 1.967  
 13473-620 Vila Bertini  
 Americana, São Paulo - Brasil  
 Tel: (+55) 19 3468-9600  
 Fax: (+55) 19 3468-9601  
 Email: [info@br.gbo.com](mailto:info@br.gbo.com)  
 Emergency phone number: (+55) 19 3468-9600

### THAILAND


Greiner Bio-One (Thailand) Ltd.  
 700/172 Moo. 1 Amata Nakorn Industrial Estate  
 Tambon Bankao, Amphur Phanthong  
 Chonburi 20160  
 Tel: (+66) 38-465633  
 Fax: (+66) 38-465635  
 Email: [office@th.gbo.com](mailto:office@th.gbo.com)  
 Emergency phone number: (+66) 38-465633

- **Recommended use / restrictions of use :**  
 To collect, transport and process blood for testing serum, plasma or whole blood in the clinical laboratory. To be used only by trained healthcare professionals according to instructions of use.

## SECTION 2: HAZARDS IDENTIFICATION

- **Classification according to NFPA 704 (Possible Rating 0-4):**

**Citric Acid Monohydrate:**  
**Health Rating:** 2  
**Flammability Rating:** 1  
**Reactivity Rating:** 0



- **Classification according to Regulation (EC) No. 1272/2008 [EU-GHS/CLP]:**

**Citric Acid Monohydrate:**  
**Eye irritation:** Category 2

- **Classification according to EU Directives 67/548/EEC or 1999/45/EC:**

**Citric Acid Monohydrate:**  
 Irritating to eyes, respiratory system and skin



- **Signal Word:**

**Citric Acid Monohydrate:** Warning



- **Hazard Statements:**

**Citric Acid Monohydrate:**  
**H319:** Causes serious eye irritation

 Page 2 / 4	<b>SAFETY DATA SHEET</b> acc. to ISO 11014, 29 CFR 1910.1200	 <b>greiner bio-one</b> A AN SOP 04.03.02-042 Rev.01 Valid from: Sept 12, 2016
	<b>VACUETTE® 9NC Coagulation Sodium Citrate 3,2%  Blood Collection Tube</b>	

- **Precautionary Statements:**  
**Citric Acid Monohydrate:**  
**P305 + P351 + P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

- **Chemical name:** Citric Acid Monohydrate  
Citric Acid Trisodium Dihydrate
- **CAS No.** Citric Acid Monohydrate 5949-29-1  
Citric Acid Trisodium Dihydrate 6132-04-03
- **Quantity of substances:** Citric Acid Monohydrate < 0,5 %  
Citric Acid Trisodium Dihydrate < 3 % } Preparation < 20 % (mass %)  
Because of trade secrets, not all components and their percentages are listed

### SECTION 4: FIRST AID MEASURES

- **Acute exposure effect:**  
May cause eye irritation. May cause mucosal irritation.
- **Skin:** Wash with soap and copious amounts of water.
- **Eyes:** Flush eyes with copious amounts of water for at least 5 minutes.
- **Inhalation:** N/A
- **Ingestion:** Dilute with large amounts of water. Consult a physician.

### SECTION 5: FIRE-FIGHTING MEASURES



- **Suitable extinguishing media:** Water spray; dry chemical; carbon dioxide.
- **Unusual fire and explosion hazards:** N/A
- **Protective equipment:** Firefighters should wear proper protective equipment and self-contained breathing apparatus with full-face piece operated in positive pressure mode.
- **Special hazards arising from the substance or mixture:** Oxides of nitrogen, carbon monoxide, carbon dioxide.

### SECTION 6: ACCIDENTAL RELEASE MEASURES

- **Personal precautions:** Avoid overexposure. Wear suitable protective clothing.
- **Methods for cleaning up:** Flush area with water.
- **Methods of containment:** Dispose as of hazardous waste. Keep in suitable, closed containers for disposal
- **Emergency procedures:** Wear respiratory protection. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas.

### SECTION 7: HANDLING AND STORAGE

- **Handling:**  
**Advice for safe handling:** Suitable for any general chemical storage area.  
**Information about protection against explosions and fires:** N/A
- **Storage:**  
Suitable for any general chemical storage area.

 greiner bio-one Page 3 / 4	<b>SAFETY DATA SHEET</b> acc. to ISO 11014, 29 CFR 1910.1200	 greiner bio-one A AN SOP 04.03.02-042 Rev.01 Valid from: Sept 12, 2016
	<b>VACUETTE® 9NC Coagulation Sodium Citrate 3,2%          Blood Collection Tube</b>	

## SECTION 8: EXPOSURE CONTROLS / PERSONAL PROTECTION



- **Control parameters:** N/A  
**Exposure limits:**  
 OSHA PEL: N/A  
 ACGIH TLV : N/A  
 Other recommended limits: N/A
- **Personal protective equipment:**  
**General protective and hygienic measures:** Avoid contact with eyes, skin and clothing.  
**Breathing equipment:** N/A  
**Hand protection:** Wear proper gloves.  
**Eye protection:** Use safety goggles.  
**Body protection:** Wear appropriate protective clothing to prevent skin exposure.
- **Hygiene measures:** Wash hands thoroughly after handling.
- **Appropriate engineering controls:** Ensure adequate ventilation. Avoid contact with skin, eyes and clothing.

## SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

- **General information**  
**Form:** Liquid (20°C)  
**Color:** Colorless  
**Odor:** Odorless
- **Odor Threshold:** N/A
- **PH-value:** 5 – 6 (20°C)
- **Change in condition**  
 Melting point/Melting range: Not determined  
 Boiling point/Boiling range: approx. 100°C
- **Flash point:** N/A
- **Flammability (solid, gaseous):** N/A
- **Danger of explosion:** Product does not present an explosion hazard
- **Vapor pressure:** Not determined
- **Density:** 1,01 g/cm<sup>3</sup>
- **Solubility in/Miscibility w/H<sub>2</sub>O:** Miscible
- **Organic solvents:** N/A
- **Solids content:** N/A
- **Partition coefficient: n-octanol/water:** N/A
- **Auto-ignition temperature:** N/A
- **Decomposition temperature:** N/A
- **Viscosity:** N/A
- **Oxidizing properties:** N/A

## SECTION 10: STABILITY AND REACTIVITY

- **Thermal decomposition / conditions to avoid:** No decomposition if used according to specifications
- **Materials to avoid:** None identified
- **Dangerous reactions:** None identified
- **Hazardous decomposition products:** Oxides of nitrogen, carbon monoxide, carbon dioxide.

 greiner bio-one Page 4 / 4	<b>SAFETY DATA SHEET</b> acc. to ISO 11014, 29 CFR 1910.1200	 greiner bio-one A AN SOP 04.03.02-042 Rev.01 Valid from: Sept 12, 2016
	<b>VACUETTE® 9NC Coagulation Sodium Citrate 3,2%          Blood Collection Tube</b>	

#### SECTION 11: TOXICOLOGICAL INFORMATION

- **Acute toxicity:** Not established
- **Eye:** May cause irritation
- **Skin:** May cause irritation
- **Inhalation:** Not determined
- **Ingestion:** None identified
- **Primary irritant effect:**
  - **On the skin:** May cause irritation
  - **On the eye:** May cause irritation
- **Sensitization:** Not established
- **Additional toxicological information:** None identified

#### SECTION 12: ECOLOGICAL INFORMATION

- **Ecotoxicological effects:** No data is available on the adverse effects of this material on the environment.
- **Other information:** N/A
- **General notes:** N/A

#### SECTION 13: DISPOSAL CONSIDERATION

- **Product: Recommendation**  
 Disposal should be done in accordance with all federal, state and local environmental regulations. Disposal must be made according to the regulations found in 40 CFR 261. This product is not a hazardous waste according to local regulations.
- **Packaging: Recommendation**  
 Disposal should be done in accordance with all federal, state and local environmental regulations. This product is not a hazardous waste according to local regulations.  
**Recommended cleansing agent:** Water if necessary with cleansing agents

#### SECTION 14: TRANSPORT INFORMATION

- **DOT regulations:** Not regulated
- **Land transport ADR/RID:** Not regulated
- **Maritime transport IMDG:** Not regulated
- **Air transport ICAO-TI and IATA-DGR:** Not regulated

#### SECTION 15: REGULATORY INFORMATION

- **OSHA, 29 CFR 1910.1200(g) and Appendix D. United Nations Globally Harmonized System of Classification and Labelling of Chemicals (GHS), third revised edition, United Nations, 2009**
- **EC directives: 91/155/EEC, 93/112/EC, 88/379/EEC**
- **Water hazard class:** N/A
- **Note:** Please note that there may be additional legal provisions to be observed. We recommend that you keep yourself informed about all applicable international, national and local regulations.

#### SECTION 16: OTHER INFORMATION

To the best of our knowledge, the information contained herein is accurate. However, neither Greiner Bio-One nor any of its subcontractors or suppliers assumes any liability whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.