

TECHNICAL DATA SHEET

Page 125 General Catalogue 2014

Polyethylene caps for Ø 15mm tubes

code 3072 and series

code	to fit tubes	bag quantity	bag weight	bag volume
3072	300705, 300900, 400705, 400900, 401100, 300707	2.500	3.20	0.016



SPECIFIC USE AND DESCRIPTION

Internal diameter: 15.4 mm
External diameter: 17.0 mm
Total length: 21.6 mm

Those three flange plug caps present significant advantages over other flange caps existing in the market, which are:

- A larger size, providing a larger surface for finger rest.
- Flanges have a larger size than other caps. More space between each flange increases the cap sealing.
- A small internal cavity designed to prevent possible drop spillage.
- A raised top for easy-on, easy-off use and to avoid cap deformation from repeated removal.
- A ribbed surface.

Note on flange cap dimensions: using flange caps implies reducing tube capacity. If as much tube volume as possible is required, we recommend the use of those caps including a chamber, such as codes 3082, 3050 and 3081. To avoid possible liquid spillage on the benchtop, leave caps stand with flanges up.

Colors to add at the end of the code:

- 02 : natural
- 05: blue
- 07: green
- 06: red

WRITTEN	VERIFIED	AUTHORISED
 David Campos Marketing Department	 Anna Mir Technical, Development and Quality Director	 Mònica Torras Sales & Marketing Director