

 DELTALAB	TECHNICAL DATA SHEET	300500 en	Date: 21/06/07
		Revision 8	Page: 1/2

Codes: 300500, 400500
Round bottom tube 11X55 mm

2007 catalogue page: 18, 20

Code	Material	Case quantity	Case weight	Case volume	Cases per pallet
300500	Crystal polystyrene	2 x 2000	7.5	0.046	40
400500	Ultraclear polypropylene	2 x 2000	5.8	0.046	40

Dimensions (+/-0.09 mm) and physical properties:


Codes:	300500	400500
Material:	Polystyrene	Polypropylene
Internal diameter (d) :	10.0 mm	9.9 mm
External diameter (D) :	11.5 mm	11.4 mm
Total length (L) :	54.4 mm	53.9 mm
Maximum volume:	3.5 ml	3.4 ml
centrifugation:	7.500 G	15.000 G



Sanitary Product. IVDD 98 / 79 / CE

SPECIFIC USE AND DESCRIPTION:

Dimensions correspond to drawing below.



Elaborated	Verified	Authorised
		

 DELTALAB	TECHNICAL DATA SHEET	300500 en	Date: 21/06/07
		Revision 8	Page: 2/2

Tubes suitable for centrifugation in normal laboratory use.

Ideal as primary tubes in applications where a small volume of blood is required, for instance in **paediatrics** and **RIA**.

Manufactured from NEW GENERATION HIGH TECHNOLOGY moulds to ensure tubes reproducibility.

One piece tubes with no joints, to ensure uniformity and dimensional accuracy of tube walls.

Caps recommended for those tubes are of the **3053** series (see description in the general catalogue).

There are many RACKS available for the tubes mentioned above. Please see our stainless steel racks code F-200 and series, Z-100, Z-500 or plastic racks code M-015, etc... Please ask for more information.

Elaborated	Verified	Authorised
		