

 deltalab	TECHNICAL DATA SHEET	19923en	Date: 25/08/10
		Revision 0	Page 1 of 2

Codes: 19923, 900025, 979930, 19924, 900028
Slide mailers




Catalogue page 2007: 160

Model	Code	Capacity	Case quantity	Weight	Volume
1	19923	1	50	0.40	0.009
2	900025	2	100	1.40	0.009
3	979930	2	250	3.50	0.040
4	19924	3	10	0.28	0.001
5	900028	5	100	1.28	0.012

SPECIFIC USE

Designed for slide transport and protection. Made of natural, autoclavable (121 °C) polypropylene. Compartment shape ensures tight fit: no risk of damaging slides. Easy snaplock. Not indicated for staining.

Elaborated	Verified	Authorised
		

	TECHNICAL DATA SHEET	19923en	Date: 25/08/10
		Revision 0	Page 2 of 2

DESCRIPTION

- Model 1 Code 19923 Container for 1 slide**

Box with a base and a hinged lid. Snap-on lid ensures a tight seal. It has a large area to write or to label. Inside the compartment, two bars of 1 mm prevent the slide from sticking to the bottom.

The container has 2 clefts that allow opening it just with one hand.

Dimensions: 80 x 40 x 6 mm.

- Model 2 Code 900025 Container for 2 slides**

Vertical container with hinged lid. The lid embodies an inside lip to prevent the slides from movement during transport. It also incorporates a hole that allows the container to be sealed in order to ensure the inviolability of the content or also to put labels.

The package is designed for 2 slides; however it can contain 4 slides.

Dimensions: 87 x 47 x 16 mm (capped).

- Model 3 Code 979930 Container for 2 slides**
Model 4 Code 19924 Container for 3 slides

Base with two or three compartments (depending on the model) and a hinged lid.

They incorporate 2 clefts on the base and the lid that allow to easy open them just with one hand.

Inside each compartment, two bars of 1 mm prevent the slide from sticking to the bottom.

Lock guides in each compartment prevent the slide from moving.

Two openings in each compartment allow catching the slide better.

Dimensions: **Code: 979930:** 84 x 71 x 6 mm (capped).

Code: 19924: 84 x 99 x 6 mm (capped).

- Model 5 Code 900028 Container for 5 slides**

Rectangular container with a hinged lid that ensures a tight seal.

Dimensions: 82 x 17 x 29 mm (capped).

See our slide line on our general catalogue.

Elaborated	Verified	Authorised
		