

	TECHNICAL DATA SHEET	440150 en	Date: 06/07/10
		Revision 3	Page: 1/1

Codes: 440150, 440151 Cytology brushes

2007 catalogue page: 148

CYTOB
CYTOBS

Code	Description	Case quantity	Case weight	Case volume
440150	Non sterile	50 x 50	4.55	0.010
440151	Sterile peel-pack	2 x 500	3.90	0.029

Instructions for use are included in sterile models.

Sanitary Product. MDD 93/42/CE.

Code 440150: Class I non sterile. 93/42/CE

Code 440151: Class I sterile. 93/42/CE

SPECIFIC USE AND DESCRIPTION

Single use product sterile by ethylene oxide.

Designed for the collection of endocervical samples for subsequent laboratory analyses (cytology, anatomy-pathology, microbiology, etc.).

The shape and material of the tip are especially designed to allow a gentle sample taking without damaging mucous membranes. The brush also allows to take enough cell samples for several subsequent clinical studies if necessary.

Product designed to facilitate the collection of endocervical samples by specialised healthcare staff and in appropriate hygienic conditions. Do not use in case of pregnancy. Do not use if packaging is broken or damaged.

Once used, product must be decontaminated or treated as a contaminated clinical residue in order to avoid cross infections. No special storage conditions are required except that integrity of packaging must be preserved.

Blue ABS shaft.

Dimensions:

Total length: 195 ± 1.5 mm

Shaft length: 175 ± 1.0 mm

Shaft diameter: 2.9 ± 0.2 mm

Tip diameter on the upper part: 5 ± 0.5 mm

Tip diameter on the lower part: 7 ± 0.5 mm

Elaborated	Verified	Approved
		