



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	162C ED
GENERAL DESCRIPTION	TAMPONE SECCO A/AL 1000PKG
DETAILED DESCRIPTION	10 boxes containing 100 tubes with aluminum stick with polyester tip and orange cap in bulk

PRODUCT COMPONENTS

APPLICATOR	ALUMINIUM APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STICK LENGTH: 150.0 mm \pm 1.0 mm TOTAL LENGTH: 152.0 mm \pm 1.0 mm STICK DIAMETER: 1.00 mm \pm 0.10 mm TIP DIAMETER: 3.00 mm \pm 0.50 mm TIP LENGTH: 15.00 mm \pm 2.00 mm
CAP	CAP YELLOW	DIMENSIONAL SPECIFICATIONS
		CAP LENGTH: 32.15 mm \pm 0.20 mm INSERTION DIAMETER: 9.30 mm \pm 0.05 mm EXTERNAL DIAMETER: 12.2 mm \pm 0.10 mm
TRANSPORT TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 155.00 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.5 mm EXTERNAL TUBE DIAMETER: 12.0 mm \pm 0.5 mm

This bulletin (Rev.00-05/09/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.

DB PR01.2



TECHNICAL DATA SHEET

PACKAGING SPECIFICATION

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	BULK

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.00-05/09/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.