



**FOR
GOOD
PIPETTING**

| | | | |
|---------|--------------------------|------------------------------|--------------------------|
| STERILE | FREE OF detectable DNase | FREE OF detectable human DNA | FREE OF detectable RNase |
| | non-cytotoxic | non-pyrogenic | |

CELLSTAR® SEROLOGICAL PIPETTES


1 ml to 50 ml

With CELLSTAR®, Greiner Bio-One offers a comprehensive range of different serological pipettes with sophisticated features. Our use of high-quality polystyrene ensures maximum transparency, with the standardised colour code allowing the volume to be quickly identified. The additional negative graduation provides an increased maximum volume with a safe distance to the filter, which protects the pipetting device against liquid aspiration. Each of our individually packed pipettes allows full traceability with the shelf life and lot number printed on the packaging.

- / For drop-free pipetting
- / Maximum accuracy
- / Clear, easy to read graduation
- / Lateral coloured Schellbach stripe on the 1, 2, 5, 10 and 25 ml pipettes simplifies volume reading

KEY FACTS

- / Sterile pipettes
- / 1 ml, 2 ml, 5 ml, 10 ml, 25 ml and 50 ml
- / 2 ml aspiration pipette without filter and graduation
- / Available in bulk or individually with plastic/plastic or paper/plastic packaging
- / Also available as triple-packed version for simplified sterile material transfer



PRODUCT OVERVIEW

FOR FURTHER INFORMATION AND/OR SAMPLE ORDERING
PLEASE VISIT OUR WEBSITE OR CONTACT US.

Serological Pipettes

Raw material: PS, Sterile: +

| Item no. | Description | Graduation | Nominal volume | Packaging | Type of packaging | Sterile | Qty. inner / outer |
|----------|---------------|------------|----------------|-------------------|-------------------|---------|--------------------|
| 604107 | 1 ml pipette | 1/100 | 1 ml | | bulk | + | 25 / 1,000 |
| 604160 | 1 ml pipette | 1/100 | 1 ml | plastic / plastic | single-packed | + | 100 / 1,000 |
| 604181 | 1 ml pipette | 1/100 | 1 ml | paper / plastic | single-packed | + | 100 / 1,000 |
| 710107 | 2 ml pipette | 1/100 | 2 ml | | bulk | + | 25 / 1,000 |
| 710160 | 2 ml pipette | 1/100 | 2 ml | plastic / plastic | single-packed | + | 100 / 1,000 |
| 710180 | 2 ml pipette | 1/100 | 2 ml | paper / plastic | single-packed | + | 100 / 1,000 |
| 606107 | 5 ml pipette | 1/10 | 5 ml | | bulk | + | 25 / 500 |
| 606160 | 5 ml pipette | 1/10 | 5 ml | plastic / plastic | single-packed | + | 50 / 200 |
| 606180 | 5 ml pipette | 1/10 | 5 ml | paper / plastic | single-packed | + | 50 / 200 |
| 607107 | 10 ml pipette | 1/10 | 10 ml | | bulk | + | 25 / 500 |
| 607160 | 10 ml pipette | 1/10 | 10 ml | plastic / plastic | single-packed | + | 50 / 200 |
| 607180 | 10 ml pipette | 1/10 | 10 ml | paper / plastic | single-packed | + | 50 / 200 |
| 760107 | 25 ml pipette | 2/10 | 25 ml | | bulk | + | 25 / 200 |
| 760160 | 25 ml pipette | 2/10 | 25 ml | plastic / plastic | single-packed | + | 50 / 200 |
| 760180 | 25 ml pipette | 2/10 | 25 ml | paper / plastic | single-packed | + | 50 / 200 |
| 768160 | 50 ml pipette | 1/2 | 50 ml | plastic / plastic | single-packed | + | 20 / 100 |
| 768180 | 50 ml pipette | 1/2 | 50 ml | paper / plastic | single-packed | + | 20 / 100 |

Serological Pipettes

Triple-packed

Raw material: PS, Packaging: plastic / plastic, Triple-packed: yes, Sterile: +

| Item no. | Description | Graduation | Nominal volume | Sterile | Qty. inner / outer |
|------------|---------------|------------|----------------|---------|--------------------|
| 604160-TRI | 1 ml pipette | 1/100 | 1 ml | + | 10 / 100 |
| 710160-TRI | 2 ml pipette | 1/100 | 2 ml | + | 10 / 100 |
| 606160-TRI | 5 ml pipette | 1/10 | 5 ml | + | 10 / 100 |
| 607160-TRI | 10 ml pipette | 1/10 | 10 ml | + | 10 / 100 |
| 760160-TRI | 25 ml pipette | 2/10 | 25 ml | + | 10 / 100 |
| 768160-TRI | 50 ml pipette | 1/2 | 50 ml | + | 10 / 100 |

Serological Pipettes

Special Models

Description: 2 ml aspiration pipette, Feature: no plug, Graduation: -, Raw material: PS, Nominal volume: 2 ml, Packaging: paper / plastic, Sterile: +

| Item no. | Sterile | Qty. inner / outer |
|----------|---------|--------------------|
| 710183 | + | 100 / 1,000 |